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DIN EN 14254: *In-vitro* diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans  
 CLSI\* GP16 "Urinalysis Approved Guideline"  
 \*CLSI (Clinical and Laboratory Standards Institute)

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## Intended use

The urine cups with screw cap are used as a specimen collection vessel for collecting, transporting and storing urine for *in vitro* diagnostics. The products are intended for use by patients (specimen collection) under the instruction of qualified medical personnel and in a professional environment by qualified medical and laboratory personnel.


## Product description

The urine cups are available in the filling volumes 70 ml and 100 ml. The urine cups consist of a clear plastic container except for one brown variant. Depending on the version, the corresponding screw caps are included or can be purchased separately.

Urine cups for 70 ml filling volume have a scale and are available optionally with or without a paper label and sterile or non-sterile.

Urine cups for 100 ml filling volume have a scale and a text field. In addition, variants with and without a paper label, double barcode label and safety label are available. Sterile and non-sterile versions are available.

## Product overview

## Urine cup with screw cap

Type	Description
<b>A</b>	Urine cup, 70 ml, scale, with and without label, sterile and non-sterile
<b>B</b>	Urine cup, 100 ml, scale, with and without label, double barcode label and safety label, sterile and non-sterile
<b>C</b>	Urine cup, 100 ml, scale, with light protection property, non-sterile

## Safety information and warnings

NOTE FOR STERILE VERSIONS: DO NOT USE THE PRODUCT IF THE BLISTER PACKAGING OR THE SAFETY LABEL IS DAMAGED.

- General precautions: Use gloves and other general personal protective equipment to protect yourself against urine and possible exposure to pathogens transmitted by biological samples.
- Handle all biological specimens and sharp collection utensils (needles) according to the guidelines and procedures of your facility.  
In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases.  
Always follow the safety guidelines and procedures of your facility.
- Always store filled urine cups closed and upright and out of the reach of children.
- The products are intended for single use. Dispose of the products and all collection utensils in containers for hazardous biological waste.
- The products may not be used after the expiry date. The expiry date of the product is the last day of the month and year as indicated.

## Storage

Store the products at room temperature.

## Transport

The products can be transported together in an ADR (packing instruction P650) conform secondary vessel.

## Limitations

1. The first morning urine is the most suitable for urine diagnostics. There should ideally be at least 4 hours between the collection of the urine sample and the last voiding.
2. Midstream urine samples are recommended over spontaneous urine because they have less contamination by urethral and/or surrounding flora. The detection of bacteria (type and number), and possibly the result of the test strip examination, can be negatively affected, by e.g. genital and/or urethral contamination.
3. The urine sample should be collected before starting antibacterial treatment if possible.
4. The storage period and temperature for urine samples depend on the expiry date of the parameters for investigation. Without additional stabilisation, urine analysis should be carried out within one to two hours.
5. The urine sample may be cooled at 2–8 °C and this can extend the stability of the sample for some analytes. There are exceptions for light-sensitive parameters such as bilirubin and urobilinogen. To protect these and other light-sensitive parameters, a urine cup 100 ml with light protection property, Art. no. 75.562.011, or a Urine Monovette® with light protection property, Art. no. 10.252.030, can be used.
6. Cooling at 2–8 °C can promote the precipitation of amorphous urates or phosphates in the urine sample and impair urine microscopy.
7. Do not freeze urine samples if they are to be used for urine sediment.
8. For cultural microbiological diagnostics, if the sample cannot be transported or immediately processed, it is recommended to store the urine either cooled at 2–8 °C or in a Urine Monovette® with boric acid, Art. no. 10.253.
9. The suitability of the urine sample for the analysis of trace elements/metals has not been tested.
10. Urine test strips are the most commonly used diagnostic aids. Before using a test strip, take note of manufacturer's information and any potential restrictions.

## Collecting and handling specimens

**READ THIS DOCUMENT COMPLETELY BEFORE STARTING URINE COLLECTION.**


### Equipment required for collecting a specimen

1. Urine cup with screw cap
2. Clean cloths or washcloth and towel

### Procedure for collecting midstream urine

The following steps are recommended for careful collection of midstream urine.

1. Clean your hands thoroughly and then clean your genital area. Dry off with a clean cloth.
2. Only for individually sterile-wrapped urine cups: Open the blister pack with the thumb and index finger using the peel-off technique. The transparent film must be facing upwards **2a**.  
Only for urine cups with a safety label (sterile): Make sure the safety label is undamaged **2b**.
3. Open the lid of the urine cup by turning it anticlockwise. Put the lid in a hygienic location. Make sure the inner side of the lid is not contaminated.
4. First void a small quantity of urine into the toilet and then, without interrupting the urine flow, start collecting the midstream urine in the cup **4b**.  
Then empty your bladder completely into the toilet **4c**.
5. Close the urine cup tightly with the screw cap by turning it clockwise. Make sure not to contaminate the inner side of the lid. Label the urine cup with your name and give the firmly closed urine cup with the sample to the responsible person immediately or proceed as instructed by your doctor.

## Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled urine cups must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

**Product-specific standards and guidelines as amended / updated**

DIN EN 14254: *In vitro* diagnostic medical devices – Single-use receptacles for the collection of specimens, other than blood, from humans  
 CLSI\* GP16 “Urinalysis Approved Guideline”  
 \*CLSI (Clinical and Laboratory Standards Institute)

**Further literature:**

European Urinalysis Guidelines: Scand J. Clin. Lab. Invest 2000; 60: 1–96.

**Key for symbols and labels:**

-  Article number
-  Batch number
-  Use by
-  CE marking
-  *In vitro* diagnostic device
-  Follow the instructions for use
-  If reused: Risk of contamination
-  Keep away from sunlight
-  Store in a dry place
-  Manufacturer
-  Country of manufacture

**The following also applies to sterile products:**

-  Do not use if package is damaged
-  Single sterile barrier system with external protective packaging
-  Sterilisation by irradiation
- or
-  Sterile fluid path
-  Do not resterilise

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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