remel Rapid SS/u System

REF R8311004..... 20 Tests/Kit

1. INTENDED USE

Remel Rapid™ SS/u System is a qualitative micromethod employing conventional and chromogenic substrates for the identification of select, medically important microorganisms commonly isolated from urine specimens. The Rapid SS/u System will provide the laboratorian with identification of common microbes associated with a positive urine culture in 2 hours. A complete listing of the organisms addressed by the Rapid SS/u System is provided in the Rapid SS/u Differential Chart.

2. SUMMARY AND EXPLANATION

The Rapid SS/u System is comprised of (1) Rapid SS/u Panels and (2) Rapid SS/u Reagent. Each Rapid SS/u Panel has several reaction cavities molded into the periphery of a plastic disposable tray. Reaction cavities contain dehydrated reactants and the tray allows the simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in Rapid Inoculation Fluid is used as the inoculum which rehydrates and initiates test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison of test results to reactivity patterns stored in the Electronic Rapid Compendium (ERIC™) database or by use of the Rapid SS/u System Differential Chart.

3. PRINCIPLE

The tests used in the Rapid SS/u System are based on microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional tests and single-substrate chromogenic tests and are described below in Table 1.

4. REAGENTS*

Rapid SS/u Reagent (provided with kit) (10 ml/Btl)
Reactive ingredient per liter:

ρ-Dimethylaminocinnamaldehyde 0.06 g

Rapid Inoculation Fluid (R8325102, supplied separately) (1 ml/Tube)

KCl 6.0 g

CaCl₂ 0.5 g

Demineralized Water 1000.0 ml

Rapid Spot Indole Reagent (R8309002, supplied separately) (15 ml/Btl)

ρ-Dimethylaminocinnamaldehyde 10.0 g

Hydrochloric Acid 100.0 ml

Demineralized Water 900.0 ml

*Adjusted as required to meet performance standards.

5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after their use. Directions should be read and followed carefully.

Caution!

- Rapid SS/u Reagent is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.
- Rapid Spot Indole Reagent may cause irritation to skin, eyes, and respiratory system.
- Refer to Safety Data Sheet for detailed information on reagent chemicals.

Composition / information on ingredients

Acetic acid 64-19-7

Hydrochloric acid 7647-01-0

2-Methoxyethanol 109-86-4

DANGER	
H315	Causes skin irritation
H319	Causes serious eye irritation
H335	May cause respiratory irritation
H336	May cause drowsiness or dizziness
H360	May damage fertility. May damage the unborn child
H373	May cause damage to organs through prolonged or repeated exposure
P201	Obtain special instructions before use
P202	Do not handle until all safety precautions have been read and understood
P281	Use personal protective equipment as required
P264	Wash face, hands and any exposed skin thoroughly after handling
P280	Wear eye/face protection
P260	Do not breathe dust/fume/gas/mist/vapors/spray
P271	Use only outdoors or in a well-ventilated area
P308+P313	IF exposed or concerned: Get medical attention/advice
P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water
P332+P313	If skin irritation occurs: Get medical advice/attention
P362	Take off contaminated clothing and wash before reuse
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P337+P313	If eye irritation persists: Get medical advice/attention
P405	Store locked up
P403+P233	Store in a well-ventilated place. Keep container tightly closed
P501	Dispose of contents/container to an approved waste disposal plant

WARNING! This product contains a chemical known in the State of California to cause birth defects or other reproductive harm.

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-

352-323-3500 (Call Collect)

6. STORAGE

Rapid SS/u System and Rapid Spot Indole Reagent should be stored in their original containers at 2-8°C until used. Allow products to equilibrate to room temperature before use. DO NOT interchange reagents among different Rapid systems. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. Rapid Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

7. PRODUCT DETERIORATION

This product should not be used if (1) the color of the reagent has changed, (2) the expiration date has passed, (3) the plastic tray is broken or the lid is compromised, or (4) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{14,15}

9. MATERIALS SUPPLIED

(1) 20 Rapid SS/u Panels, (2) 20 report forms, (3) Rapid SS/u Reagent (one plastic dropper-bottle containing reagent sufficient for 20 panels), (4) 2 chipboard incubation trays, (5) Instructions for use (IFU).

10. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Gram stain reagents, (7) Microscope slides, (8) Oxidase reagent, (9) Cotton swabs, (10) Rapid Inoculation Fluid-1 ml (R8325102), (11) McFarland #1 turbidity standard or equivalent (R20411), (12) Pipettes, (13) Rapid Spot Indole Reagent (R8309002), (14) ERIC (Electronic Rapid Compendium, R8323600).

11. CONTENTS SYMBOLS

SS/u Panels	SS/u Panels
Rapid Report Forms	Rapid Report Forms
SS/u Reagent	SS/u Reagent
Incubation Trays	Incubation Trays

12. PROCEDURE

Inoculum Preparation:

- Only urine isolates should be selected for testing.** The use of isolates from other body sites or fluids is not recommended.

Notes:

- Test isolate colonial morphology should be closely examined, since mixed microbial populations cannot be used as inoculum. Where there is an indication of polymicrobial isolation, each colony type should be isolated and processed in a Rapid SS/u panel independently.
- Where appropriate, isolates should be examined by Gram stain, wet mount, or oxidase test prior to use in the system.
- Test organisms may be removed from a variety of selective and nonselective agar growth media. The following types of media are recommended: Tryptic Soy Agar (TSA) with or without 5% Sheep Blood; Eosin Methylene Blue (EMB) Agar; Phenylethyl Alcohol (PEA) Agar; Nutrient Agar; MacConkey Agar.

Notes:

- Beta-hemolytic streptococci should not be tested using the Rapid SS/u System. The Rapid STR System is recommended for these isolates.
 - Some media containing or supplemented with mono- or disaccharides (e.g., Sabouraud Dextrose Agar or Mannitol Salt Agar) are not recommended since they may suppress glycolytic activity and reduce test selectivity.
 - Plates used for inoculum preparation should preferably be 18-24 hours old. Slow-growing isolates may be tested using 48-hour plates.
 - The use of media other than those recommended may compromise test performance.
- Using a cotton swab or inoculating loop, suspend sufficient growth from the agar plate culture in Rapid Inoculation Fluid (1 ml) to achieve a visual turbidity equal to a #1 McFarland turbidity standard or equivalent.

Notes:

- Suspensions significantly less turbid than a #1 McFarland standard will result in aberrant reactions.
 - Suspensions slightly more turbid than a #1 McFarland standard will not affect test performance and are recommended for stock cultures and quality control strains. However, suspensions prepared with a turbidity far greater than a #1 McFarland standard will compromise test performance.
 - Suspensions should be mixed thoroughly and vortexed if required.
 - Suspensions should be used within 15 minutes of preparation.
- An agar plate may be inoculated for purity and any additional testing that may be required using a loopful of the test suspension from the inoculation fluid tube. Incubate the plate for 18-24 hours at 35-37°C.

Inoculation of Rapid SS/u Panels:


- Peel back the lid of the panel over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
- Using a pipette, gently transfer the entire contents of the Inoculation Fluid tube into the upper right-hand corner

Table 1. Principles and Components of the Rapid SS/u System


Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
Before Reagent Addition:					
1	GMS	Amino acid arylamide	0.5%	Hydrolysis of the peptide substrate releases free yellow p nitrophenol.	1
2	ONPG	o-Nitrophenyl-β,D-galactoside	0.25%		
3	G1	p-Nitrophenyl-β,D-glycoside	0.25%		
4	G2	p-Nitrophenyl-β,D-glycoside	0.25%	Hydrolysis of the colorless nitrophenylated glycoside releases yellow σ- or p nitrophenol.	2-7
5	G3	p-Nitrophenyl-β,D-glycoside	0.25%		
6	PHS	p-Nitrophenyl phosphate	0.5%	Hydrolysis of the colorless phosphoester releases yellow p-nitrophenol.	2
7	URE	Urea	0.9%	Hydrolysis of urea produces basic products which raise the pH and change the indicator.	2
After Reagent Addition:					
7	IND	Tryptophane	0.5%	Utilization of tryptophane results in the formation of indole which is detected with Rapid Spot Indole Reagent.	2
8	A1	Amino acid-β-naphthylamide	0.05%	Hydrolysis of the aryl-substituted amide releases β-naphthylamine which is detected with Rapid SS/u Reagent.	1, 8-13
9	A2	Amino acid-β-naphthylamide	0.05%		
10	A3	Amino acid-β-naphthylamide	0.05%		

of the panel. Reseal the inoculation port of the panel by pressing the peel-back tab back in place.

- After adding the test suspension, and while keeping the panel on a level surface, tilt the panel back away from the reaction cavities at approximately a 45-degree angle (see below).


- While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.

- While maintaining a level, horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction cavities until the inoculum flows along the baffles into the reaction cavities (see below). This should evacuate all of the inoculum from the rear portion of the panel.

Note: If the panel is tilted too quickly, air may be trapped at the test cavity junction, restricting fluid movement.

- Return the panel to a level position. If necessary, gently tap the panel on the bench top to remove any air trapped in the cavities.

Notes:

- Examine the test cavities, which should appear bubble-free and uniformly filled. Slight irregularities in test cavity fills are acceptable and will not affect test performance. If the panel is grossly misfilled, a new panel should be inoculated and the misfilled panel discarded.
- Complete the inoculation of each panel receiving inoculation fluid before inoculating additional panels.
- Do not allow the inoculum to rest in the back portion of the panel for prolonged periods without completing the procedure.

Incubation of Rapid SS/u Panels:

Incubate inoculated panels at 35-37°C in a non-CO₂ incubator for 2 hours. For ease of handling, panels may be incubated in the chipboard incubation trays provided with the kit.

Rapid SS/u Panel Test Location

Cavity #	1	2	3	4	5	6	7	8	9	10
Test Code	GMS	ONPG	G1	G2	G3	PHS	URE	A1	A2	A3
							IND	Rapid SS/u Reagent		

Table 2. Interpretation of Rapid SS/u System Tests*

Cavity #	Test Code	Reagent	Reaction		Comments
			Positive	Negative	
Before Reagent Addition:					
1	GMS	None	Medium or bright yellow	Clear, tan, or very pale yellow	Only the development of a distinct yellow color should be scored as positive. Pale yellow or a hint of yellow should be scored as negative.
2	ONPG				
3	G1				
4	G2				
5	G3				
6	PHS				
7	URE	None	Dark red, red, red-orange, or orange	Yellow	Any development of a red or orange color should be scored as positive.
After Reagent Addition:					
7	IND	Rapid Spot Indole Reagent	Black, brown, or tan	Orange	Any development of a tan, brown, or muddy color should be scored as positive.
8	A1	Rapid SS/u Reagent	Purple, violet, red, or dark pink	Yellow, orange, or pale pink	Only significant color development should be scored as positive. Pale shades of color should be scored as negative.
9	A2				
10	A3				

*NOTE: Panels should be read by looking down through the reaction cavities against a white background.

Table 3. Quality Control Chart for RapID SS/u Panels

Organism	GMS	ONPG	G1	G2	G3	PHS	URE	A1	A2	A3	IND
<i>Escherichia coli</i> ^a ATCC® 25922	+	+	+	-	-	-	-	-	V	-	+
<i>Enterococcus faecalis</i> ^a ATCC® 29212	-	-	-	-	+	-	-	V	-	+	-
<i>Klebsiella pneumoniae</i> ATCC® 13883	+	+	-	+	-	-	V	-	+	V	-
<i>Proteus mirabilis</i> ATCC® 29906 or 25933	+	-	-	-	-	+	+	-	+	-	-
<i>Serratia marcescens</i> ATCC® 29882 or 8100	+	V	-	-	+	V	-	+	+	V	-

+, positive; -, negative; V, variable

^aKey indicator strains demonstrate acceptable performance of the most labile substrate in the system and reactivity in a significant number of wells, according to Clinical and Laboratory Standards Institute recommendations for streamlined quality control.²³

Notes:

- The quality control of RapID reagents is accomplished by obtaining the expected reactions for tests requiring the addition of the reagents (cavities 7-10).
- Organisms which have been repeatedly transferred on agar media for prolonged periods may provide aberrant results.
- Quality control strains should be stored frozen or lyophilized. Prior to use, quality control strains should be transferred 2-3 times from storage on an agar medium that is recommended for use with the RapID SS/u System.
- Formulations, additives, and ingredients of culture media vary from manufacturer to manufacturer and may vary from batch to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve quality control discrepancies.

15. LIMITATIONS

- The use of RapID SS/u System and the interpretation of results requires the knowledge of a competent laboratorian who is trained in general microbiological methods and judiciously makes use of training, experience, specimen information, and other pertinent procedures before reporting the identification obtained using this system.
- Characteristics such as Gram stain reaction, oxidase, and cellular and colonial morphology must be considered when using the RapID SS/u System.
- The RapID SS/u System must be used with pure cultures of test organisms. The use of mixed microbial

populations or direct testing of clinical material without culture will result in aberrant results.

- The RapID SS/u System is designed for use with common urine isolates including the taxa listed in the RapID SS/u Differential Chart. The use of isolates from other body sites or organisms not specifically listed in the chart may lead to misidentifications.
- Expected values listed for RapID SS/u System tests may differ from conventional test results or previously reported information.
- The accuracy of the RapID SS/u System is based upon the statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single test found in the RapID SS/u System to establish the identification of a test isolate is subject to the error inherent in that test alone.

16. PERFORMANCE CHARACTERISTICS

The RapID SS/u System performance characteristics have been established by laboratory testing of reference and stock cultures at Remel and by clinical evaluations using fresh clinical and stock isolates.^{16,17}

17. BIBLIOGRAPHY

- Cerney, G. 1978. Eur. J. Appl. Microbiol. 5:113-122.
- Blazevic, D.J. and G.M. Ederer. 1975. Principles of Biochemical Tests in Diagnostic Microbiology. John Wiley & Sons, New York, NY.
- Brisou, B., C. Richard, and A. Leuroit. 1972. Ann. Microbiol. l'Institute Pasteur (Paris). 123:341-347.
- Ewing, W.H. 1986. Edwards and Ewing's Identification of Enterobacteriaceae. 4th ed. Elsevier Science Publishing Company, Inc., New York, NY.
- Farmer III, J.J., B.R. Davis, F.W. Hickman-Brenner, A. McWhorter, G.P. Huntley-Carter, M.A. Asbury, C.

Riddle, H.G. Wathen-Grady, C. Elias, G.R. Fanning, A.G. Steigerwalt, C.M. O'Hara, G.K. Morris, P.B. Smith, and D. J. Brenner. 1985. J. Clin. Microbiol. 21:46-76.

- Halstead, D.C., M.R. Hoffert, and G.G. Colasante. 1987. J. Clin. Microbiol. 25:42-44.
- Holt J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams. 2000. Bergey's Manual of Determinative Bacteriology. 9th ed. Lippincott Williams and Wilkins., Philadelphia, PA.
- Giammanco, G., J. Bussiere, M. Toucos, G. Brault, and L. LeMinor. 1980. Ann. Microbiol. l'Institute Pasteur (Paris). 131A:181-187.
- Lenette, E.H., A. Balows, W.J. Hausler, Jr., and H.J. Shadomy. 1985. Manual of Clinical Microbiology. 4th ed. ASM, Washington, D.C.
- Nord, C.E., A.A. Lindberg, and A. Dahlback. 1975. Med. Microbiol. Immunol. 161:231-238.
- Peterson, E.H. and E.J. Hsu. 1978. J. Food Sci. 43:1853-1856.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Westley, J.R., P.J. Anderson, V.A. Close, B. Halpern, and E.M. Lederberg. 1967. Appl. Microbiol. 15-822-825.
- Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C..
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- DeGirolami, P.C., J.S. Poliferno, L.S. Mills, and K.A. Eichelberger. 1988. Am. J. Clin. Pathol. 89:791-793.
- Morganstern, F. A. and S.B. Potter. 1988. Lab. Med. 19:368-370.
- Grimont, P.A.D., F. Grimont, and H.L.C. Dulong de Rosnay. 1977. J. Gen. Microbiol. 98:39-66.
- Hansen, W. and E. Yourassowsky. 1984. J. Clin. Microbiol. 20:1177-1179.
- Kilian, M. and P. Bulow. 1976. Acta. Pathol. Microbiol. Scand. B. 84:245-251.
- Mulczyk, M. and A. Szewczuk. 1970. J. Gen. Microbiol. 61:9-13.
- Mutjens, H.L., van der Ros-van de Repe, and J. van Druten. 1984. J. Clin. Microbiol. 20:684-686.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

18. PACKAGING

REF R8311004 RapID SS/u System..... 20 Tests/Kit

19. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	For Laboratory Use Only
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Manufactured by

RapID™ is a trademark of Thermo Fisher Scientific and its subsidiaries.

ERIC™ is a trademark of Thermo Fisher Scientific and its subsidiaries.

ATCC® is a registered trademark of American Type Culture Collection.

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DANGER	H315	Causes skin irritation
	H319	Causes serious eye irritation
US ONLY	H335	May cause respiratory irritation
	H336	May cause drowsiness or dizziness
US & EU	H360	May damage fertility. May damage the unborn child
	H373	May cause damage to organs through prolonged or repeated exposure
	P201	Obtain special instructions before use
	P202	Do not handle until all safety precautions have been read and understood
	P281	Use personal protective equipment as required
	P264	Wash face, hands and any exposed skin thoroughly after handling
	P280	Wear eye/face protection
	P260	Do not breathe dust/fume/gas/mist/vapors/spray
	P271	Use only outdoors or in a well-ventilated area
	P308+P313	IF exposed or concerned: Get medical attention/advice
	P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
	P302+P352	IF ON SKIN: Wash with plenty of soap and water
	P332+P313	If skin irritation occurs: Get medical advice/attention
	P362	Take off contaminated clothing and wash before reuse
	P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
	P337+P313	If eye irritation persists: Get medical advice/attention
	P405	Store locked up
	P403+P233	Store in a well-ventilated place. Keep container tightly closed
	P501	Dispose of contents/container to an approved waste disposal plant

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RapID™ è un marchio di Remel Inc.
ERIC™ è un marchio di Remel Inc.