www.thermofisher.com

Mannitol Salt Agar

REF PO5027A

Intended Use

Mannitol Salt Agar (PO5027A) is a selective medium for the isolation of pathogenic staphylococci from clinical samples, including wound, faecal and respiratory tract samples. Mannitol Salt Agar is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having staphylococcal infections.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Staphylococci are Gram-positive cocci, widespread in nature and generally non-pathogenic^{1,2}. However, some staphylococci can cause disease in humans^{1,2}. The staphylococcus most frequently associated with disease in humans is *Staphylococcus aureus*¹, an opportunistic pathogen commonly found as a commensal of the human upper respiratory tract and skin microbiomes². *S. aureus* enters the body through cuts or abrasions where it can cause a spectrum of disease, from minor infections such as cellulitis and abscesses to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis².

Principle of Method

Mannitol Salt Agar devices contain Lab Lemco and peptone as a source of nitrogen and other nutrients. Mannitol fermentation aids the differentiation of *Staphylococcus* spp. and phenol red is the pH indicator; it is yellow at acid pH and red/pink at alkaline pH. High salt content will inhibit most other organisms. Sodium chloride at 7.5% will inhibit many non-target organisms.

Typical Formula*

	grams per litre
'Lab-Lemco' powder	1.0
Peptone	10.0
Mannitol	10.0
Sodium chloride	75.0
Phenol red	0.025
Agar	15.0

Physical Appearance

Colour	Antique pink
Clarity	Transparent
Fill weight	17 ± 5%
pH	7.5 ± 0.2

Materials Provided

10 x 90mm Mannitol Salt agar plates.

Each plate should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI).

Procedure

The method followed may depend upon the sample under test. Follow local protocols and guidelines. For clinical samples:

- Inoculate and streak the specimen onto the medium.
- If material is being cultured directly from a swab, roll the swab over a small surface area and streak for the isolation of single colonies.
- Incubate the medium aerobically at 36°C ± 1°C for 36 to 48 hours.
- Suspect *Staphylococcus aureus* colonies are yellow with yellow zones. Mannitol negative staphylococci are small colonies surrounded by pink zones.
- Inoculate suspect colonies onto a non-selective medium before identifying using appropriate methods. Testing of other sample types such as food may require a different incubation protocol.

Interpretation

The presence of yellow colonies with yellow halos indicates *Staphylococcus aureus* strains. The presence of pink

colonies with pink medium indicates *Staphylococcus epidermidis*.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 36 - 48 h @ 36 ± 1 °C ° ± 2°C

Positive Controls Inoculum level: 50-120 cfu (colony forming units) Colony count is ≥ 50% of the control medium (tryptone soya agar) count.	
<i>Staphylococcus aureus</i> ATCC® 6583™	1-2mm, yellow shiny colonies.
<i>Staphylococcus aureus</i> ATCC® 25923™	1-2mm, yellow shiny colonies.
<i>Staphylococcus epidermidis</i> ATCC® 12228™	1mm, white shiny colonies.
Negative Controls Inoculum level: >10 ⁴ cfu	
<i>Escherichia coli</i> ATCC® 25922™	No growth.

Limitations

Identifications are presumptive and colonies must be confirmed using appropriate biochemical/serological methods. It is recommended that strains are sub-cultured to a non-selective medium containing less salt before carrying out identification tests, for example the tube coagulase test. Some strains of *S. aureus* may exhibit a delayed mannitol fermentation reaction; negative plates should be re-incubated for a further 24 hours. Extended incubation may result in accompanying flora growing through. Enterococci may grow as pin-point colonies. Strains of *S. aureus* with particular growth requirements may not grow on the medium and those with atypical enzyme patterns may give anomalous results.

Performance Characteristics

Performance was evaluated using 17 bacterial strains including the following: *S. aureus*, *S. epidermidis*, *S. saprophyticus*, *Enterococcus faecalis*, *E. coli* and *Proteus mirabilis*. All organisms gave expected growth characteristics according to the current product specification.

Bibliography

- Public Health England. 2020a. "Identification of Staphylococcus Species, Micrococcus Species and Rothia Species." UK Standards for Microbiology Investigations. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
- Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.

Symbol Legend

Symbol	Definition
	Catalogue number

	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
Made in Germany	Made in Germany

ATCC Licensed Derivative[®]

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

Version	Modifications Introduced
1.0	14/07/07 - Original Document

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Oxoid Deutschland GmbH, Am Lippeglacis