www.thermofisher.com

MacConkey Agar No. 3 Modified

REF PO5140A & PB5224E*

*This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Columbia CNA Aesculin Selective Agar (product code: PB5115A)

This IFU is intended for side 2 – MacConkey Agar No 3 Modified side of the biplate PB5224E.

Intended Use

MacConkey Agar No 3 (Modified) (PO5140A) is a selective medium for the isolation of *Acinetobacter* spp. and differentiation between lactose and non-lactose fermenting Gram-negative bacteria with the inhibition of Gram-positive cocci from a wide range of clinical samples, such as urine and faeces, or from swabs taken from wounds or infections.

Columbia CNA Aesculin Selective Agar / MacConkey Agar No. 3, mod. (PB5224E) is a biplate that enables isolation of *Acinetobacter* spp. and differentiation between lactose and non-lactose fermenting Gram-negative bacteria with the inhibition of Gram-positive cocci from a wide range of clinical samples, such as urine and faeces, or from swabs taken from wounds or infections.

Columbia CNA Aesculin Selective Agar / MacConkey Agar No. 3, mod. (PB5224E) and MacConkey Agar No 3 (Modified) (PO5140A) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having lactose and non-lactose fermenting Gram-negative bacterial infections. Columbia CNA Aesculin Selective Agar / MacConkey Agar No. 3, mod. (PB5224E) and MacConkey Agar No 3 (Modified) (PO5140A) can also be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Acinetobacter* spp. infection.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Pathogenic aerobic and anaerobic microorganisms encompass a broad range of microorganisms, which can often be found as part of the normal skin and mucosal flora of humans and the infections from aerobic microorganisms can have a range of manifestations ranging from superficial to life-threatening (Public Health England; 2015). *Escherichia coli* (*E. coli*) are non-spore forming bacteria that are able to grow in aerobic and anaerobic condition (Public Health England, 2015).

Salmonella is a bacterial pathogen that can be isolated from faeces, blood, bone marrow, bile, urine, food, animal feed and environmental materials. Ingestion of contaminated food and water can cause foodborne infections, including gastroenteritis, typhoid fever, paratyphoid fever or even death in humans. All *Salmonella* serotypes can cause disease in human (Public Health England, 2021).

Acinetobacter baumannii is a Gram-negative nosocomial pathogen that can persist on dry surfaces longer than any other Gram-negative bacteria. Its ability to persist in moist and dry surfaces for more than 20 days

Thermo

SCIENTIFIC

contributes to its widespread in a hospital setting (Gerner-Smith, P. 1995). *Acinetobacter* spp. are commonly isolated from locations such as hand, groin, toe webs etc (Public Health England 2015). Due to the high antibiotic resistance shown by this bacterium an early identification is often recommended. *Acinetobacter* spp. have been isolated in connection with community-acquired and nosocomial pneumonias, urogenital tract, eye and soft tissue infections and are difficult to treat particularly due to their broad antibiotic resistance

Principle of Method

MacConkey Agar No.3 Modified is more selective when compared to the original MacConkey Agar due to the inclusion of specially prepared fraction of bile salts and crystal violet. Peptone supplies nutrients and agar is the solidifying agent. Bile salts are inhibitory to non-intestinal bacteria and help to prevent the swarming of *Proteus* spp. The bile salt mixture in the MacConkey agar no. 3 Modified was modified to allow the growth of *Acinetobacter* spp. as well. Bile salts and crystal violet inhibit the growth of Gram-positive cocci. Lactose is added as a carbon source. Differentiation of bacteria is achieved by the combination of lactose and the indicator dye neutral red, which is red at acid pH and yellow at alkaline. Bacteria which ferment lactose appear as red - pink-coloured colonies which may be surrounded by zones of precipitated bile salts. Precipitation is caused by the action of the acid produced by lactose fermentation on the bile salts. Bacteria which do not ferment lactose, such as *Salmonella*, usually appear as colourless to straw colonies. Therefore, MacConkey agar No. 3 Modified can be used with clinical specimens likely to contain mixed flora such as urine and wounds as it allows the preliminary grouping of Gram-negative bacteria into lactose fermenters and non-fermenters.

Typical Formula

	grams per litre
Peptone	20.0
Lactose	10.0
Bile salts No. 3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar	15.0

Physical Appearance

Colour	Antique pink
Clarity	Transparent
Total Fill weight	17 ± 5%
pH	7.1 ± 0.2

Materials Provided

PO5140A: 10 x 90mm MacConkey agar No.3 Modified plates

PB5224E: 10 x 90mm Columbia CNA Aesculin Selective Agar with Blood / MacConkey Agar No. 3 Modified biplates

Each plate should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

The performance of this medium can be verified by testing the following reference strains.

Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 1, B 4, B 11, B 14, B 17, B 25, B26, B37, B57 and B60.

Procedure

- Allow product to equilibrate to room temperature.
- Inoculate and streak the specimen onto the medium using a standard loop.
- **PO5140A:** Incubate plates aerobically for 18–48 hours at 36 ± 1°C.
- **PB5224E:** Incubate plates aerobically for 24 hours at 36 ± 1°C.
- Visually inspect plates to assess colony growth and colour under good lighting.

Interpretation

The presence of pink colonies and bile precipitation indicates *Escherichia coli*. The presence of straw colonies indicates *Salmonella* Typhimurium. Rose shiny colonies on (PO5140A) or Brownish shiny colonies on (PB5224E) indicates *Acinetobacter baumannii*

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

PO5140A:

Incubation Conditions: 18 - 48 h @ 36 ± 1°C, aerobic

Positive Controls Inoculum level: 50-120 cfu Colony count is ≥ 50% of the control medium count.	
<i>Escherichia coli</i> ATCC® 8739™	3 – 7 mm, pink colonies, bile precipitate
<i>Salmonella</i> Typhimurium ATCC® 14028™	3 - 4 mm, brown shiny colonies
<i>Acinetobacter baumannii</i> ATCC® 19606™	Good growth, brownish shiny colonies (40h)
Inoculum level 10 ³ -10 ⁴ cfu	
<i>Proteus mirabilis</i> ATCC® 29906™	Good growth, brownish shiny colonies
Negative Controls Inoculum level: ≥10 ⁴ cfu	
<i>Staphylococcus aureus</i> ATCC® 6538	Partial inhibition (≤ 100 cfu)

PB5224E:

Incubation Conditions: 24 h @ 36 + 1°C aerobic

Positive Controls Inoculum level: 10 ³ -10 ⁴ cfu	
<i>Acinetobacter baumannii</i> ATCC® 19606™	Brownish, shiny colonies
<i>Escherichia coli</i> ATCC® 25922™	Pink colonies and bile precipitation
Negative Controls Inoculum level: 10 ⁴ -10 ⁵ cfu	
<i>Enterococcus faecalis</i> ATCC® 29212™	No growth

Limitations

Organisms with atypical enzyme patterns may give anomalous reactions on MacConkey Agar No 3 (Modified).

All identifications are presumptive and should be confirmed using appropriate methods. Prolonged incubation may lead to confusing results; therefore, do not incubate beyond 48 hours. Some strains with specific growth requirements may grow poorly or not at all on this medium. Organisms which are resistant to the selective agents may be able to grow. The swarming of proteus is inhibited on this medium but occasional strains may be able to swarm.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of lactose and non-lactose fermenting Gram negative bacteria and *Acinetobacter* spp. is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of Columbia CNA Aesculin Selective Agar / MacConkey Agar No. 3, mod. (PB5224E) was demonstrated by an overall pass rate of 100% obtained for the product over 2 months of testing (20.04.2022 – 03.06.2022). The precision of MacConkey Agar No 3 (Modified) (PO5140A) was demonstrated by an

overall pass rate of 100% obtained for the product over twelve months of testing (28.04.2021- 19.04.2022). This shows that the performance is reproducible.

MacConkey Agar No 3 (Modified) (PO5140A) device is tested in-house as part of the QC process since the products were launched in 2003. For target organisms, when using 50-120 cfu inoculum of *Escherichia coli*, *Salmonella* Typhimurium and *Acinetobacter baumannii* and incubating the device at: 18 – 48 h @ 36 ± 1°C, aerobic, conditions., the user can recover organisms with colony size and morphology as listed in this document. For target organism *Proteus mirabilis* using an inoculum level of 10³ – 10⁴ cfu and incubating at 18 – 48 h @ 36 ± 1°C, aerobic conditions., the user can recover organisms with colony size and morphology that meets the defined acceptance criteria For non-target organisms, when using ≥ 10⁴ cfu of *Staphylococcus aureus* and incubating the device at : 18 – 48 h @ 36 ± 1°C, aerobic, the user can expect partial inhibition or no growth.

Columbia CNA Aesculin Selective Agar / . MacConkey Agar No. 3, mod. (PB5224E) device is tested in-house as part of the QC process since the products were launched in 2003. For target organisms, when using 10³ – 10⁴ cfu inoculum of *Acinetobacter baumannii* and *Escherichia coli* and incubating the device at: 24 h @ 36 ± 1°C, aerobic, conditions., the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10⁴ – 10⁵ cfu inoculum of *Enterococcus faecalis* and incubating the device at 24 h @ 36 ± 1°C, aerobic, the user can expect no growth.

Bibliography

1. Gerner-Smidt, P. (1995). Taxonomy and epidemiology of *Acinetobacter* infections. Rev Med Microbiol., 6, 186-195.
2. Public Health England. 2015. "Identification of Enterobacteriaceae." UK Standards for Microbiology Investigations. UK SMI ID 16 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
3. Public Health England. 2021. "Identification of Salmonella species." UK Standards for Microbiology Investigations UK SMI ID 24 Issue 4.
<https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit

	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in Germany

ATCC Licensed Derivative[®]

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

Version	Date of modifications introduced
1.0	2022.09.14. Original Document

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Bibliografia

1. Gerner-Smidt, P. (1995). Taxonomy and epidemiology of *Acinetobacter* infections. *Rev Med Microbiol.*, 6, 186-195.
2. Public Health England. 2015. "Identification of Enterobacteriaceae." UK Standards for Microbiology Investigations. UK SMI ID 16 Issue 4. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
3. Public Health England. 2021. "Identification of Salmonella species." UK Standards for Microbiology Investigations UK SMI ID 24 Issue 4.

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Oxoid Deutschland GmbH, Am Lippeglacis