Mueller Hinton Agar with Sheep Blood

REF PB5007A

Intended Use

Mueller Hinton Agar with Sheep Blood is recommended for disc diffusion susceptibility testing as standardized by the Clinical and Laboratory Standards Institute (CLSI) to aid clinicians in testing fastidious organisms isolated from clinical samples.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Un-supplemented Mueller-Hinton Agar (MHA) is suitable for testing rapidly growing aerobic bacteria when using disc or gradient strip diffusion methods; however, it is not adequate when testing more fastidious microorganisms. MHA with 5% defibrinated sheep blood is recommended by the Clinical and Laboratory Standards Institute (CLSI) for testing of bacteria such as streptococci and *Neisseria meningitidis*^{1,2,3}. Mueller-Hinton Agar was designed to be a reproducible culture medium for the isolation of pathogenic *Neisseria* species⁴. Specific media for *Neisseria* spp. have since replaced Mueller-Hinton Agar and its major use is now antimicrobial susceptibility testing.

Principle of Method

Mueller-Hinton media contain beef extract, casein hydrolysate to provide nitrogen, carbon, vitamins, amino acids, sulphur and other essential nutrients. Starch is added to absorb toxic metabolites that are produced and its hydrolysis yields dextrose, which serves as an energy source. Agar is added as the solidifying agent.

Mueller-Hinton Agar is a stable, reproducible, and standardized solid medium which is recommended for the study of susceptibility of bacteria to antimicrobial agents, using the Kirby-Bauer method of disc diffusion or dilution in agar.

Unsupplemented Mueller-Hinton Agar is recommended by CLSI and EUCAST for antimicrobial susceptibility testing of non-fastidious strains. It is also recommended in the Deutsches Institut für Normung (DIN) standard.

Supplements (such as defibrinated blood and nicotinamide adenine dinucleotide (NAD)) are added to make Mueller-Hinton Agar suitable for more fastidious strains.

Typical Formula

	grams per litre
Beef, dehydrated infusion from	300.0
Casein hydrolysate	17.5
Starch	1.5
Agar	17.0

Additions

Defibrinated sheep blood	50.0ml
--------------------------	--------

Physical Appearance

Colour	Signal red
Clarity	Opaque
Fill weight	23.5g ± 5%
pH	7.4 ± 0.2

Materials Provided

PB5007A: 10 x 90mm Mueller Hinton Agar with Sheep Blood plates

Each plate should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q 5.

Procedure

Follow the procedures in the appropriate standard. After incubation, measure the zone of inhibition as recommended and interpret the result in accordance with the method being followed for example CLSI^{1,2,3}.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium,

and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

This medium can be tested with the following strain in accordance with CLSI methodology

Incubation Conditions: 20 - 24 h @ 36 ± 1°C in an enhanced carbon dioxide atmosphere

Target Organism	Zone size range (mm)
<i>Streptococcus pyogenes</i> ATCC® 49619™	
Moxifloxacin 5	25-31
Erythromycin 15	25-30
Cefpodoxime 10	28-34
Ertapenem 10	28-35

Limitations

The medium may not support the growth of strains with specific growth requirements. Close adherence to recognised standard methodology is required to obtain satisfactory results.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of fastidious microorganisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Mueller Hinton Agar with Sheep Blood (PB5007A) was demonstrated by an overall pass rate of 100% obtained for the product over two months of testing (23.JUN.2022–16.AUG.2022; 10 batches). This shows that the performance is reproducible.


Mueller Hinton Agar with Sheep Blood (PB5007A) is tested in-house as part of the QC process since 2003. The inoculum is prepared from working cultures, diluted to obtain the specified inoculum levels and the test media is inoculated. Antibiotic discs are loaded in the same order and concentrations as they appear in the inspection plan for the medium. Following incubation at 36 ± 1°C for 20-24 hours in an enhanced carbon dioxide atmosphere, the zone of inhibition diameter of the test organisms is measured, and the results and conclusions can be created.

Bibliography

- (1) CLSI: Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; 3rd Ed CLSI Guideline M45 Wayne PA, Clinical and Laboratories Standards Institute 2016
- (2) CLSI: Performance Standards for Antimicrobial Disk Susceptibility Tests 13th Ed CLSI standard M02 Wayne PA, Clinical and Laboratories Standards Institute 2018
- (3) CLSI: Performance standards for antimicrobial susceptibility testing. 30th Ed CLSI supplement M100 Wayne PA, Clinical and Laboratories Standards Institute 2020
- (4) Mueller J. H. and Hinton J (1941) Proc. Soc. Exp. Biol. and Med. 48. 330-333.

Symbol Legend

Symbol	Definition
	Catalogue number

	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in Germany

ATCC Licensed Derivative®

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. CLSI is a trademark of the Clinical Laboratory and Standards Institute. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

Version	Date of modifications introduced
1.0	2022-09-29 Original Document

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Zdroje

(1) CLSI: Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; 3rd Ed CLSI Guideline M45 Wayne PA, Clinical and Laboratories Standards Institute 2016

(2) CLSI: Performance Standards for Antimicrobial Disk Susceptibility Tests 13th Ed CLSI standard M02 Wayne PA, Clinical and Laboratories Standards Institute 2018

(3) CLSI: Performance standards for antimicrobial susceptibility testing. 30th Ed CLSI supplement M100 Wayne PA, Clinical and Laboratories Standards Institute 2020

(4) Mueller J. H. and Hinton J (1941) Proc. Soc. Exp. Biol. and Med. 48. 330-333.

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©2022 Thermo Fisher Scientific Inc. Todos los derechos reservados.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus filiales.
Oxid Deutschland GmbH, Am Lippeglacis