www.thermofisher.com

Thermo
SCIENTIFIC

Lactose	10.0
L-Cysteine	0.128
Bromothymol blue	0.02
Agar	15.0

CLED Agar

REF PO5009A & PO5217E*

*This instructions for use (IFU) document is intended to be read in conjunction with the IFU for MacConkey Agar No.3 (product code: PO5002A)

This IFU is intended to be used for Side 1 (CLED Agar) of the CLED Agar/ MacConkey Agar No.3 biplate.

Intended Use

CLED Agar and CLED Agar / MacConkey Agar No. 3 (PO5009A and PO5217E) devices are differential media for the isolation and enumeration of the common micro-organisms causing urinary tract infections from urine samples.

CLED Agar (PO5009A and PO5217E) devices are used with urinary samples in a diagnostic workflow for pre-natal screening and to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated, nor a companion diagnostic.

Summary and Explanation

Urinary tract infections, often referred to as UTIs, are common infections that occur when bacteria invade the urinary tract by entering via the urethra. These infections can impact different parts of the urinary tract, but the most prevalent is cystitis. UTIs occur more frequently in women, as the urethra is shorter, and therefore closer to the rectum. This increases the likelihood of infection as it is easier for bacteria to enter the urinary tract. Kidney infections are another form of UTI which is called acute pyelonephritis (APN) These are more serious than cystitis, but far less common.

The symptoms of cystitis can include frequent urination, or needing to urinate whilst having an empty bladder, a burning sensation or pain whilst urinating, blood in the urine and pressure or cramps in the lower abdominal area. The symptoms of kidney infections are very different. These can include chills, a fever, nausea and vomiting and pain in the lower back¹. There are many species of bacteria that cause UTIs, but the three most common are *Escherichia coli*, *Staphylococcus aureus* and *Proteus mirabilis*.

Principle of Method

Peptones and Lab Lemco powder are present to supply the compounds and peptides required for the growth of bacteria and agar is the solidifying agent. Lactose provides a carbohydrate source. Bromothymol blue is a pH indicator which differentiates lactose fermenters (yellow) from non-fermenters. Cystine enhances the growth of cystine-dependent coliforms. Electrolytes are reduced in order to restrict the swarming of *Proteus* spp.

Typical Formula

	<u>grams per litre</u>
Peptone	4.0
'Lab-Lemco' powder	3.0
Tryptone	4.0

Physical Appearance

Colour	Pale Green
Clarity	Transparent
Fill weight	17g ± 5%
pH	7.3 ± 0.2

Materials Provided

PO5009A: 10 x 90mm CLED agar plates
PO5217E: 10 x 90mm CLED agar/ MacConkey Agar No.3 agar biplates

Each plate should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 41 and Q 5.

Agar/ MacConkey Agar No.3 (Side 1) PO5217E and incubating the device at $36 \pm 1^\circ\text{C}$ for 18-24 hours, the user can recover organisms with colony size and morphology as listed in this document

Procedure

- Allow product to equilibrate to room temperature.
- Inoculate and streak the specimen onto the medium using a standard loop.
- Incubate plates aerobically for 18–24 hours at $36 \pm 1^\circ\text{C}$.
- Visually inspect plates to assess colony growth and colour under good lighting.

Interpretation

- Light shiny yellow colonies indicate *Escherichia coli*
- Yellow shiny colonies indicate *Staphylococcus aureus*
- Light green shiny colonies indicate *Proteus vulgaris*

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 - 24 h @ $36 \pm 1^\circ\text{C}$ aerobic

Positive Controls	
Inoculum level 50 – 120 cfu Colony count is $\geq 50\%$ of the control medium count.	
<i>Escherichia coli</i> ATCC® 25922™	2 - 4mm, light yellow shiny colonies
<i>Staphylococcus aureus</i> ATCC® 25923™	2 – 3mm, yellow shiny colonies
Inoculum level $10^3 - 10^4$ cfu	
<i>Proteus vulgaris</i> ATCC® 8427™	Good growth, light green shiny colonies

Limitations

Organisms with atypical enzyme patterns may give anomalous reactions on CLED Agar.

This medium is for laboratory use only and is not intended for the diagnosis of disease or other conditions. Identifications are presumptive and colonies should be identified using appropriate methods. The medium is non-selective but some strains with specific growth requirements may not grow for example *Shigella* species may not grow on electrolyte-deficient media.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of UTI- causing organisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of CLED Agar (PO5009A) and CLED Agar / MacConkey Agar No. 3 (PO5217E) was demonstrated by an overall pass rate of 100% obtained for the product over 14 years of testing 2008-2022; 654 batches for PO5009A and 575 batches for PO5217E. This shows that the performance is reproducible.

For target organisms, when using an inoculum level of 50-120 cfu of *Staphylococcus aureus* or *Escherichia coli* and $10^3 - 10^4$ cfu of *Proteus vulgaris* on CLED Agar (PO5009A) and incubating the device at $36 \pm 1^\circ\text{C}$ for 18-24 hours, the user can recover organisms with colony size and morphology as listed in this document

When using an inoculum level of $10^3 - 10^4$ cfu of *Staphylococcus aureus*, and *Escherichia coli* on CLED

Bibliography

1. Centers for Disease Control and Prevention. 2021. 'Urinary Tract Infection'. <https://www.cdc.gov/antibiotic-use/uti.html>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier

ATCC Licensed Derivative®

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

Version	Modifications Introduced
1.0	2022-09-01 Original Document

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Oxoid Deutschland GmbH, Am Lippeglacis 4-8,
©2022 Thermo Fisher Scientific Inc. Toate drepturile rezervate.
Fisher Scientific Inc. și a filialelor sale.