Oxoid Prepared Medium

Oxoid C.L.E.D. Medium / MacConkey Agar No. 3

REF PO5217E

Intended Use

IVD

Thermo Scientific™ Oxoid™ C.L.E.D. Medium / MacConkey Agar No. 3 (biplate) is used for the microbiological examination of urine samples.

For professional use only (in vitro diagnostic use).

Summary and Explanation

C.L.E.D. Medium

Cystine Lactose Electrolyte Deficient (C.L.E.D.) Medium is made to the formula described by Mackey and Sandys¹ as a modification for urinary bacteriology of the Electrolyte Deficient Medium developed by Sandys². This medium is recommended for urinary bacteriology, supporting the growth of all urinary pathogens and giving good colonial differentiation and clear diagnostic characteristics. The presence of important contaminants such as diphtheroids, lactobacilli and micrococci are also clearly elicited, giving an indication of the degree of contamination. In the laboratory C.L.E.D. Medium provides a valuable non-inhibitory diagnostic agar for plate culture of urinary organisms. It is electrolyte deficient to prevent the swarming of *Proteus* species. C.L.E.D. Medium may also be used for wound swabs.

MacConkey Agar No. 3

MacConkey agar No. 3 is based on the original formulation of MacConkey^{1,2}. It is used for isolating and differentiating lactose fermenting from lactose non-fermenting Gram-negative enteric bacilli. This is a more selective modification of MacConkey agar which is suitable for the detection and enumeration of coliform organisms and also for the detection and isolation of *Salmonella* and *Shigella* species occurring in pathological and food specimens. Due to the inclusion of a specially prepared fraction of bile salts in addition to crystal violet, the medium gives improved differentiation between coliforms and lactose non-fermenting organisms whilst gram-positive cocci are inhibited. This formulation corresponds with that recommended by the American Public Health Association³ for the direct plating of water samples for coliforms, for the examination of food samples for food poisoning organisms⁴ and for the isolation of *Shigella* species in food⁵. Amongst other examples of the use of MacConkey Agar no 3 are the enumeration of coli-aerogenes in animal faeces⁶, poultry faecal and carcass samples^{7,8} and canned minced chicken⁹. MacConkey agar No. 3 can be used with clinical specimens likely to contain mixed flora such as urine and wounds as it allows the preliminary grouping of Gram-negative bacteria into lactose fermenters and non-fermenters¹⁰.

Principle

C.L.E.D. Medium

Peptones and Lab Lemco powder are present to supply the compounds and peptides required for the growth of bacteria and agar is the solidifying agent. Lactose provides a carbohydrate source. Bromothymol blue is a pH indicator which differentiates lactose fermenters (yellow) from non-fermenters. Cystine enhances the growth of cystine-dependent coliforms. Electrolytes are reduced in order to restrict the swarming of *Proteus* spp.

Typical Colony Appearance (18 hours incubation)

<i>Escherichia coli</i>	Yellow, opaque, slightly deeper coloured centre. Non lactose fermenting strains appear as blue
<i>Klebsiella</i> and <i>Enterobacter</i>	Yellow to blue mucoid
<i>Salmonella</i>	Blue flat
<i>Enterococcus faecalis</i>	Yellow, approximately

	0.5 mm diameter
<i>Staphylococcus aureus</i>	Yellow colonies approximately 0.75 mm
Coagulase negative staphylococci	Pale yellow to white
Corynebacteria	Very small grey colonies
Lactobacilli	Similar to corynebacteria but with a rougher surface

MacConkey Agar No. 3

Peptone supplies nutrients and agar is the solidifying agent. Bile salts are inhibitory to non-intestinal bacteria and help to prevent proteus swarming. Bile salts and crystal violet inhibit the growth of Gram-positive cocci. Lactose is added as a carbon source. Differentiation of bacteria is achieved by the combination of lactose and the indicator dye neutral red which is red at acid pH and yellow at alkaline. Bacteria which ferment lactose appear as red-pink coloured colonies which may be surrounded by zones of precipitated bile salts. Precipitation is caused by the action of the acid produced by lactose fermentation on the bile salts. Bacteria which do not ferment lactose, such as salmonellae, usually appear as colourless to straw colonies.

Typical Formula*

C.L.E.D. Medium:

	grams per litre
Peptone	4.0
'Lab Lemco' powder	3.0
Tryptone	4.0
Lactose	10.0
L-Cystine	0.128
Bromothymol blue	0.02
Agar	15.0

Side 2 - MacConkey Agar No. 3: grams per litre

Peptone	20.0
Lactose	10.0
Bile salts No.3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar	15.0

* Adjusted as required to meet performance standards

Physical Characteristics

Side 1 - C.L.E.D. Medium:

Colour	Pale green
Clarity	Transparent
pH	7.3 ± 0.2

Side 2 - MacConkey Agar No. 3:

Colour	Antique pink
Clarity	Transparent
pH	7.1 ± 0.2

Total fill weight 17g ± 5%

Precautions

This product is for *in vitro* diagnostic use and should only be used by trained individuals. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully.

For professional use only. Safety Data Sheet available on request.

Storage

This product is ready to use and no further preparation is necessary. Store product in its original packaging at 2–12°C until used.

Store away from light.
Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following the recommended guidelines.

Materials Required but Not Supplied

- (1) Inoculating loops, swabs, collection containers
- (2) Incubators
- (3) Quality control organisms

More information can be found at www.thermofisher.com

Procedure

Refer to any local guidelines and protocols. The method used will depend upon the specimen type.

For urine samples:

- (1) Inoculate the specimen as soon as possible after it is received in the laboratory. The use of a calibrated loop is recommended to facilitate counting colonies and estimate the colony forming units (cfu) /mL.
- (2) Streak the sample onto the C.L.E.D. medium first and then onto the MacConkey No. 3 agar.
- (3) Incubate aerobically at 36 ± 1°C for 18 to 24 hours.
- (4) Determine colony count and observe colonial characteristics and morphology.

Quality Control

This medium can be tested with the following strains:

Incubation Conditions: 18 – 24 h @ 36 ± 1°C, aerobic

Side 1 - C.L.E.D. Medium:

Positive Controls	
Inoculum of 10 ³ – 10 ⁴ colony forming units (cfu).	
<i>Escherichia coli</i> ATCC®25922	Good growth, yellow, opaque colonies.
<i>Staphylococcus aureus</i> ATCC®25923	Good growth, deep yellow colonies.

Side 2 - MacConkey Agar No. 3:

Positive Controls	
Inoculum of 10 ³ – 10 ⁴ colony forming units (cfu).	
<i>Escherichia coli</i> ATCC®25922	Good growth, pink colonies and bile precipitate.
<i>Salmonella</i> Typhimurium ATCC®14028	Good growth, brownish irregular colonies.
Negative Control	
Inoculum of ≥10 ⁴ cfu.	
<i>Enterococcus faecalis</i> ATCC®29212	No growth.

Note:

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The product should not be used if:

- (1) The product is contaminated
- (2) The colour has changed

- (3) The expiration date has passed
- (4) There are other signs of deterioration

Performance

Performance was evaluated using 29 bacterial strains including the following; *E.coli*, *Klebsiella* spp, *Proteus* spp, *Salmonella* Typhimurium, *Salmonella* Nottingham, *S. aureus*, *S. epidermidis*, *S. saprophyticus*, *Enterococcus* spp. and *Corynebacterium* spp. All organisms gave expected growth characteristics according to the current product specification.

Limitations

C.L.E.D. Medium

Identifications are presumptive and colonies should be identified using appropriate methods. The medium is non-selective but some strains with specific growth requirements may not grow, for example *Shigella* species may not grow on electrolyte-deficient media.

MacConkey Agar No. 3

Some strains with specific growth requirements may grow poorly or not at all on this medium. Organisms which are resistant to the selective agents may be able to grow. The swarming of proteus is inhibited on this medium but occasional strains may be able to swarm.

Packaging

PO5217E	Ten 90mm bi-plates, film wrapped
---------	----------------------------------

Bibliography

C.L.E.D. Medium

- (1) Mackey J. P. and Sandys G. H. (1966) *B.M.J.* 1. 1173.
- (2) Sandys G. H. (1960) *J. Med. Lab. Techn.* 17. 224.

MacConkey Agar No. 3

- (1) MacConkey A.T., 1900 *The Lancet*, Volume 156, (4010) 20
- (2) MacConkey A.T., 1905 *J Hyg*, Jul; 5(3): 333–379.
- (3) American Public Health Association (1998) *Standard Methods for the Examination of Water and Wastewater. 20th Edn. APHA Inc. Washington DC.*
- (4) American Public Health Association (1976) *Compendium of methods for the Microbiological Examination of Foods. APHA Inc. Washington DC.*
- (5) ISO 21567:2004 *Microbiology of food and animal feeding stuffs - Horizontal method for the detection of Shigella spp.*
- (6) Medrek T. F. and Barnes E. M. (1962) *J. Appl. Bact.* 25(2)159-168.
- (7) Barnes E. M. and Goldberg H. S. (1962) *J. Appl. Bact.* 25(1) 94-106.
- (8) Barnes E.M. and Shrimpton D.H. (1957) *J. Appl. Bact* 20(2) 273-285
- (9) Thornley M.J. (1957) *J. Appl. Bact.* 20(2) 286-298.
- (10) Forbes, B.A., Sahm D. F. and Weissfeld A. S. 1998. *Bailey & Scott's diagnostic microbiology*, 10th ed. Mosby-Year Book, Inc., St. Louis.

Symbol Legend

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation (storage temp.)
	Use by (expiration date)

	Lot number
	Protect from light
	Consult instructions for use

 The ATCC Licensed Derivative[®] Emblem, the ATCC Licensed Derivative[®] word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC[®] cultures.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. ATCC[®] is a trademark of American Type Culture Collection. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

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
MacConkey Agar No. 3

- (1) MacConkey A.T., 1900 *The Lancet*, Volume 156, (4010) 20
- (2) MacConkey A.T., 1905 *J Hyg*, Jul; 5(3): 333-379.
- (3) American Public Health Association (1998) *Standard Methods for the Examination of Water and Wastewater. 20th Edn. APHA Inc. Washington DC.*
- (4) American Public Health Association (1976) *Compendium of methods for the Microbiological Examination of Foods. APHA Inc. Washington DC.*
- (5) ISO 21567:2004 *Microbiology of food and animal feeding stuffs - Horizontal method for the detection of Shigella spp.*
- (6) Medrek T. F. and Barnes E. M. (1962) *J. Appl. Bact.* 25(2)159-168.
- (7) Barnes E. M. and Goldberg H. S. (1962) *J. Appl. Bact.* 25(1) 94-106.
- (8) Barnes E.M. and Shrimpton D.H. (1957) *J. Appl. Bact.* 20(2) 273-285
- (9) Thornley M.J. (1957) *J. Appl. Bact.* 20(2) 286-298.
- (10) Forbes, B.A., Sahm D. F. and Weissfeld A. S. 1998. *Bailey & Scott's diagnostic microbiology*, 10th ed. Mosby-Year Book, Inc., St. Louis.

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Agar MacConkey N.º 3

- (1) MacConkey A.T., 1900 *The Lancet*, Volume 156, (4010) 20
- (2) MacConkey A.T., 1905 *J Hyg*, Jul; 5(3): 333-379.
- (3) American Public Health Association (1998) *Standard Methods for the Examination of Water and Wastewater. 20th Edn. APHA Inc. Washington DC.*
- (4) American Public Health Association (1976) *Compendium of methods for the Microbiological Examination of Foods. APHA Inc. Washington DC.*
- (5) ISO 21567:2004 *Microbiology of food and animal feeding stuffs - Horizontal method for the detection of Shigella spp.*
- (6) Medrek T. F. and Barnes E. M. (1962) *J. Appl. Bact.* 25(2)159-168.
- (7) Barnes E. M. and Goldberg H. S. (1962) *J. Appl. Bact.* 25(1) 94-106.

- (8) Barnes E.M. and Shrimpton D.H. (1957) *J. Appl. Bact* 20(2) 273-285
 (9) Thornley M.J. (1957) *J. Appl. Bact.* 20(2) 286-298.
 (10) Forbes, B.A., Sahn D. F. and Weissfeld A. S. 1998. *Bailey & Scott's diagnostic microbiology*, 10th ed. Mosby-Year Book, Inc., St. Louis.

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Gélose MacConkey N° 3 :

- (1) MacConkey A. T., 1900 *The Lancet*, Volume 156, (4010) 20
- (2) MacConkey A.T., 1905 *J Hyg*, Jul; 5(3): 333-379.
- (3) American Public Health Association (1998) *Standard Methods for the Examination of Water and Wastewater. 20th Edn. APHA Inc. Washington DC.*
- (4) American Public Health Association (1976) *Compendium of methods for the Microbiological Examination of Foods. APHA Inc. Washington DC.*
- (5) ISO 21567:2004 *Microbiology of food and animal feeding stuffs - Horizontal method for the detection of Shigella spp.*
- (6) Medrek T. F. and Barnes E. M. (1962) *J. Appl. Bact.* 25(2)159-168.
- (7) Barnes E. M. and Goldberg H. S. (1962) *J. Appl. Bact.* 25(1) 94-106.
- (8) Barnes E.M. and Shrimpton D.H. (1957) *J. Appl. Bact* 20(2) 273-285
- (9) Thornley M. J. (1957) *J. Appl. Bact.* 20(2) 286-298.

(10) Forbes, B.A., Sahn D. F. and Weissfeld A. S. 1998. *Bailey & Scott's diagnostic microbiology*, 10th ed. Mosby-Year Book, Inc., St. Louis.

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
MacConkey Agar No. 3

- (1) MacConkey A.T., 1900 *The Lancet*, Volume 156, (4010) 20
- (2) MacConkey A.T., 1905 *J Hyg*, Jul; 5(3): 333-379.
- (3) American Public Health Association (1998) *Standard Methods for the Examination of Water and Wastewater. 20th Edn. APHA Inc. Washington DC.*
- (4) American Public Health Association (1976) *Compendium of methods for the Microbiological Examination of Foods. APHA Inc. Washington DC.*
- (5) ISO 21567:2004 *Microbiology of food and animal feeding stuffs - Horizontal method for the detection of Shigella spp.*
- (6) Medrek T. F. and Barnes E. M. (1962) *J. Appl. Bact.* 25(2)159-168.
- (7) Barnes E. M. and Goldberg H. S. (1962) *J. Appl. Bact.* 25(1) 94-106.
- (8) Barnes E.M. and Shrimpton D.H. (1957) *J. Appl. Bact.* 20(2) 273-285
- (9) Thornley M.J. (1957) *J. Appl. Bact.* 20(2) 286-298.

(10) Forbes, B.A., Sahn D. F. and Weissfeld A. S. 1998. *Bailey & Scott's diagnostic microbiology*, 10th ed. Mosby-Year Book, Inc., St. Louis.

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Versione 3

Meio preparado Oxoid

Oxoid C.L.E.D. Medium/MacConkey Agar No. 3

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MacConkey Agar No. 3

- (1) MacConkey A.T., 1900 *The Lancet*, Volume 156, (4010) 20
- (2) MacConkey A.T., 1905 *J Hyg*, Jul; 5(3): 333-379.
- (3) American Public Health Association (1998) *Standard Methods for the Examination of Water and Wastewater. 20th Edn. APHA Inc. Washington DC.*
- (4) American Public Health Association (1976) *Compendium of methods for the Microbiological Examination of Foods. APHA Inc. Washington DC.*
- (5) ISO 21567:2004 *Microbiology of food and animal feeding stuffs - Horizontal method for the detection of Shigella spp.*
- (6) Medrek T. F. and Barnes E. M. (1962) *J. Appl. Bact.* 25(2)159-168.
- (7) Barnes E. M. and Goldberg H. S. (1962) *J. Appl. Bact.* 25(1) 94-106.
- (8) Barnes E.M. and Shrimpton D.H. (1957) *J. Appl. Bact* 20(2) 273-285
- (9) Thornley M.J. (1957) *J. Appl. Bact.* 20(2) 286-298.

(10) Forbes, B.A., Sahm D. F. and Weissfeld A. S. 1998. *Bailey & Scott's diagnostic microbiology*, 10th ed. Mosby-Year Book, Inc., St. Louis.

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Oxoid Deutschland GmbH