MacConkey Agar No. 3

**REF PO5002A & PB5207E*,
PO5217E****

*This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Columbia Agar with Blood (product code: PB5008A)

**This instructions for use (IFU) document is intended to be read in conjunction with the IFU for CLED Agar (product code: PO5009A)

This IFU is intended for Side 2 - MacConkey Agar No.3 side of the biplates PB5207E and PO5217E.

Intended Use

MacConkey Agar No. 3 (PO5002A) device is a selective medium for the isolation and differentiation between lactose and non-lactose fermenting Gram-negative bacteria with the inhibition of Gram-positive cocci from a wide range of clinical samples, such as urine and faeces, or from swabs taken from wounds or infections.

Columbia Agar with Blood / MacConkey Agar No. 3 (PB5207E) is a biplate that enables the isolation and differentiation of lactose fermenting and non-lactose fermenting Gram-negative bacteria with the inhibition of Gram-positive cocci from a wide range of clinical samples, such as urine and faeces, or from swabs taken from wounds or infections.

CLED Agar / MacConkey Agar No. 3 (PO5217E) is a biplate for the isolation and differentiation between lactose and non-lactose fermenting Gram-negative bacteria with the inhibition of Gram-positive cocci from a wide range of clinical samples, such as urine and faeces, or from swabs taken from wounds or infections.

Columbia Agar with Blood / MacConkey Agar No. 3 (PB5207E), MacConkey Agar No. 3 (PO5002A) & CLED Agar / MacConkey Agar No. 3 (PO5217E) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having lactose and non-lactose fermenting Gram-negative bacterial infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Pathogenic aerobic and anaerobic microorganisms encompass a broad range of microorganisms, which can often be found as part of the normal skin and mucosal flora of humans and the infections from aerobic microorganisms can have a range of manifestations ranging from superficial to life-threatening¹. *Escherichia coli* are non-spore forming bacteria that are able to grow in aerobic and anaerobic conditions¹.

Salmonella is a bacterial pathogen that can be isolated from faeces, blood, bone marrow, bile, urine, food, animal feed and environmental materials. Ingestion of contaminated food and water can cause foodborne infections, including gastroenteritis, typhoid fever, paratyphoid fever or even death in humans. All *Salmonella* serotypes can cause disease in human².

Principle of Method

MacConkey agar No. 3 is based on the original formulation of MacConkey Agar. MacConkey Agar No.3 is more selective when compared to the original MacConkey Agar due to the inclusion of specially prepared fraction of bile salts and crystal violet. Peptone supplies nutrients and agar is the solidifying agent. Bile salts are inhibitory to non-intestinal bacteria and help to prevent the swarming of *Proteus* spp. Bile salts and crystal violet inhibit the growth of Gram-positive cocci. Lactose is added as a carbon source. Differentiation of bacteria is achieved by the combination of lactose and the indicator dye neutral red, which is red at acid pH and yellow at alkaline. Bacteria which ferment lactose appear as red - pink-coloured colonies which may be surrounded by zones of precipitated bile salts. Precipitation is caused by the action of the acid produced by lactose fermentation on the bile salts. Bacteria which do not ferment lactose, such as *Salmonella*, usually appear as colourless to straw colonies. Therefore, MacConkey agar No. 3 can be used with clinical specimens likely to contain mixed flora such as urine and wounds as it allows the preliminary grouping of Gram-negative bacteria into lactose fermenters and non-fermenters.

Typical Formula

	grams per litre
Peptone	20.0
Lactose	10.0
Bile salts No. 3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar	15.0

Physical Appearance

PO5002A:

Colour	Antique pink to heather violet
Clarity	Clear
Fill weight	19 ± 0.2g
pH	7.1 ± 0.2

PB5207E & PO5217E:

Colour	Antique pink
Clarity	Transparent
Total fill weight	17 ± 5%
pH	7.1 ± 0.2

Materials Provided

PO5002A: 10 x 90mm MacConkey agar No.3 plates
PB5207E: 10 x 90mm Columbia Agar with Blood / MacConkey Agar No. 3 biplates
PO5217E: 10 x 90mm CLED Agar / MacConkey Agar No. 3 biplates

Each plate should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store products PB5207E and PO5217E in their original packaging at 2–12°C until used.
- Store products PO5002A in the original packaging at 2–25°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.

- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 1, B 4, B 11, B 14, B 17, B 25, B26, B37, B57 and B60.

Procedure

- Allow product to equilibrate to room temperature.
- Inoculate and streak the specimen onto the medium using a standard loop.
- **PO5002A:** Incubate plates aerobically for 18–24 hours at 30 ± 35 °C.
- **PB5207E & PO5217E:** Incubate plates aerobically for 18–24 hours at 36 ± 1°C.
- Visually inspect plates to assess colony growth and colour under good lighting.

Interpretation

The presence of pink colonies and bile precipitation indicates *Escherichia coli*. The presence of straw colonies indicates *Salmonella* Typhimurium.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

PO5002A:

Incubation Conditions: 18 - 24 h @ 30 - 35°C aerobic

Positive Controls

Inoculum level: 50-120 cfu Colony count is ≥ 50% of the control medium count.	
<i>Escherichia coli</i> ATCC® 8739	Pink colonies and bile precipitation
<i>Salmonella</i> Typhimurium ATCC® 14028	Straw colonies
Negative Controls Inoculum level: 10 ⁴ -10 ⁵ cfu Incubation Conditions: 18 - 48 h @ 32 + 1°C aerobic	
<i>Staphylococcus aureus</i> ATCC® 6538	No Growth

PB5207E & PO5217E:

Incubation Conditions: 18 - 24 h @ 36 + 1°C aerobic

Positive Controls Inoculum level: 10 ³ -10 ⁴ cfu	
<i>Escherichia coli</i> ATCC® 25922™	Good growth, pink colonies and bile precipitation
<i>Salmonella</i> Typhimurium ATCC® 14028™	Good growth, straw colonies
Negative Controls Inoculum level: ≥10 ⁴ cfu	
<i>Enterococcus faecalis</i> ATCC® 29212™	No growth

Limitations

Organisms with atypical enzyme patterns may give anomalous reactions on MacConkey Agar No 3.

All identifications are presumptive and should be confirmed using appropriate methods. Prolonged incubation may lead to confusing results; therefore, do not incubate beyond 48 hours. Some strains with specific growth requirements may grow poorly or not at all on this medium. Organisms which are resistant to the selective agents may be able to grow. The swarming of proteus is inhibited on this medium but occasional strains may be able to swarm.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of lactose and non-lactose fermenting Gram negative bacteria and *Acinetobacter* spp. is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of Columbia Agar with Blood / MacConkey Agar No. 3 (PB5207E) was demonstrated by an overall pass rate of 100% obtained for the product over 2 months of testing (20.04.2022 – 02.06.2022; 10 batches). The precision of CLED Agar / MacConkey Agar No. 3 (PO5217E). was demonstrated by an overall pass rate of 100% obtained for the product over 2 months of testing (07.03.2022 – 23.05.2022). The precision of MacConkey Agar No. 3 (PO5002A). was demonstrated by an overall pass rate of 100% obtained for the product over 2 months of testing (19.05.2022-02.06.2022). This shows that the performance is reproducible.

MacConkey Agar No. 3 (PO5002A) device is tested in-house as part of the QC process since the products were launched in 2003. For target organisms, when using 50-120 cfu inoculum of *Escherichia coli* and *Salmonella* Typhimurium, and incubating the device at: 18 – 24 h @ 30 – 35°C, aerobic, conditions., the user can recover

organisms with colony size and morphology as listed in this document. For non-target organisms, when using $10^4 - 10^5$ cfu *Staphylococcus aureus* and incubating at $18 - 48$ h @ $32 \pm 1^\circ\text{C}$, aerobic, conditions., the user can expect no growth.

Columbia Agar with Blood / MacConkey Agar No. 3(PB5207E) device is tested in-house as part of the QC process since the products were launched in 2003. For target organisms, when using 10^3-10^4 cfu inoculum of *Escherichia coli* and *Salmonella* Typhimurium and incubating the device at $18 - 24$ h @ $36 \pm 1^\circ\text{C}$, aerobic conditions., the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using $\geq 10^4$ cfu inoculum of *Enterococcus faecalis* and incubating the device at $18 - 24$ h @ $36 \pm 1^\circ\text{C}$, aerobic the user can expect no growth.


CLED Agar / MacConkey Agar No. 3 (PO5217E) device is tested in-house as part of the QC process since the products were launched in 2003. For target organisms, when using $10^3 - 10^4$ cfu inoculum of *Escherichia coli* and *Salmonella* Typhimurium ATCC® 14028™ and incubating the device at: $18 - 24$ h @ $36 \pm 1^\circ\text{C}$, aerobic, conditions., the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using $\geq 10^4$ cfu of *Enterococcus faecalis* and incubating the device at $18 - 24$ h @ $36 \pm 1^\circ\text{C}$, aerobic, the user can expect no growth.

Bibliography

- Public Health England. 2015a. " Identification of Enterobacteriaceae." UK Standards for Microbiology Investigations. UK SMI ID 16 Issue 4.. <https://www.gov.uk/government/publications/smi-id-16-identification-of-enterobacteriaceae>
- Public Health England. 2021a. "Introduction to the preliminary identification of medically important bacteria and fungi from culture." UK Standards for Microbiology Investigations. UK SMI ID 1 Issue 3 Accessed 19th January 2022. <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use

	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in UK
	Made in Germany

ATCC Licensed Derivative[®]

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

Version	Date of modifications introduced
1.0	2022.09.14. Original Document

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Oxoid Deutschland GmbH, Am Lippeglacis