www.thermofisher.com

Thermo
SCIENTIFIC

Ox bile 20.0
Agar 10.0

Enterococcus Selective Agar (Bile Aesculin Azide Agar)

REF PO5062A

Intended Use

Enterococcus Selective Agar (Bile aesculin Azide Agar) is a selective medium for the isolation and presumptive identification of enterococci from clinical samples (e.g. wounds, throat, sputum, faeces, etc.). Enterococcus Selective Agar is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Enterococci are Gram-positive microorganisms found in the gastrointestinal tract of humans and animals, plants, soil, water, environment and fermented products¹. The most common species isolated from human infections are *E. faecalis* and *E. faecium*¹. Enterococcus species are significant human pathogens and therefore are under surveillance by the European Antimicrobial Resistance Surveillance Network (EARS-Net)². Enterococci demonstrate resistance to a number of antibiotics including β -lactams, aminoglycosides, chloramphenicol and last-line antibiotics such as vancomycin¹. The percentage of *E. faecium* vancomycin resistance between 2014 to 2020 has increased dramatically from 9% to 14%². It has been estimated that 1.27 million deaths were attributable to antimicrobial resistance in pathogens such as *E. faecalis* and *E. faecium* in 2019. Additionally, the number of reported invasive *Enterococcus faecium* isolates, mostly bloodstream infections, increased at EU/EEA level in 2020 compared to 2019².

Principle of Method

Enterococcus Selective Agar is used for the isolation, presumptive identification and enumeration of enterococci. Enterococci, including *E. faecalis* and *E. faecium*, were previously classified as Group D streptococci and are still commonly referred to as Group D streptococci. These bacteria can hydrolyse aesculin into aesculetin and glucose. The reaction between aesculetin and ferric ammonium citrate results in a brown-black to black halo around the colonies. Ox bile inhibits the growth of Gram-positive bacteria except enterococci, while sodium azide suppresses the Gram-negative bacteria.

Tryptone in the media provides a nitrogen source, as well as other nutrients to support the growth of the target bacteria. Yeast extract supplies other growth factors, such as vitamins. Enterococcus Selective Agar contains Aesculin Azide Agar (CM0591B) as a base medium.

Typical Formula

| | grams per litre |
|-------------------------|-----------------|
| Tryptone | 20.0 |
| Yeast extract | 5.0 |
| Sodium chloride | 5.0 |
| Aesculin | 1.0 |
| Ferric ammonium citrate | 0.5 |
| Sodium azide | 0.55 |
| Sodium citrate | 1.0 |

Physical Appearance

Colour Sand yellow
Clarity Transparent
Fill weight 17.0 g \pm 5%
pH 7.0 \pm 0.2

Materials Provided

PO5062A: 10 x 90mm Enterococcus Selective Agar (Bile Aesculin Azide Agar) plates

Each plate should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and applicable local regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following recommended local guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 4 and Q 5.

Procedure

- Allow product to equilibrate to room temperature.

2. Organisation for Economic Co-operation and Development. 2022. 'Antimicrobial Resistance in the EU/EEA – A One Health response'. <https://www.ecdc.europa.eu/en/publications-data/antimicrobial-resistance-eueea-one-health-response>

- Inoculate and streak the specimen onto the medium using a standard loop.
- Incubate plates aerobically for 18–24 hours at 36 ± 1 °C.
- Visually inspect plates to assess colony growth and colour under good lighting.

Interpretation

Brown colonies with brown halos indicates *Enterococcus faecalis*. Brown shiny colonies with brown halos indicates *Enterococcus faecium*.

Quality Control

It is the responsibility of the user to perform Quality Control testing, taking into account the intended use of the medium, and in accordance with any applicable local regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 - 24 h @ 36 ± 1°C aerobic

| | |
|--|--|
| Positive Controls Inoculum level: 50 – 120 cfu Colony count is ≥ 50% of the control medium count. | |
| Strain tested by membrane filtration | |
| <i>Enterococcus faecalis</i> ATCC® 29212™ | Brown colonies and brown halos |
| Strain tested by spread plate | |
| <i>Enterococcus faecium</i> ATCC® 19434™ | 1-2 mm, brown shiny colonies with brown halos. |
| Negative Controls Inoculum level: ≥ 10 ⁴ cfu | |
| <i>Escherichia coli</i> ATCC® 25922™ | Total inhibition |
| <i>Streptococcus pyogenes</i> ATCC® 12344™ | Total inhibition |

Limitations

Strains of the target organisms with particular growth requirements or which are sensitive to the selective agents in the medium may grow weakly or not at all. Further tests may be required to confirm the identity of the isolates.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of Enterococci species is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of Enterococcus Selective Agar (Bile Aesculin Azide Agar) (PO5062A) was demonstrated by an overall pass rate of 100% obtained for the product over ten months of testing (18-Aug-2021 – 14-Jun-2022; ten recent batches). This shows that the performance is reproducible.

Enterococcus Selective Agar (Bile Aesculin Azide Agar) (PO5062A) is tested in-house as part of the QC process since 2003. For target organisms, when using 50-120 cfu inoculum of *Enterococcus faecalis* or *Enterococcus faecium* by incubating the device at 36 ± 1°C for 18-24 hours the user can recover organisms as listed in this document.

Bibliography

1. Public Health England. 2021. 'Identification of *Streptococcus* species, *Enterococcus* species

Symbol Legend

| Symbol | Definition |
|--------|---|
| | Catalogue number |
| | <i>In Vitro</i> Diagnostic Medical Device |
| | Batch code |
| | Temperature limit |
| | Use-by date |
| | Keep away from sunlight |
| | Do not re-use |
| | Consult instructions for use or consult electronic instructions for use |
| | Contains sufficient for <n> tests |
| | Do not use if packaging damaged and consult instructions for use |
| | Manufacturer |
| | Authorized representative in the European Community/ European Union |
| | European Conformity Assessment |
| | UK Conformity Assessment |
| | Unique device identifier |

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of
American Type Culture Collection.
All other trademarks are the property of Thermo Fisher
Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

| Version | Date of modifications introduced |
|---------|----------------------------------|
| 1.0 | 2022-11-03 Original Document |

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Deutschland GmbH, Am Lippeglacis