Oxoid Prepared Medium

Oxoid Columbia Agar with Sheep Blood^{PLUS}

REF PB5039A

Intended Use

IVD

Thermo Scientific™ Oxoid™ Columbia Agar with Sheep Blood^{PLUS} is a multi-purpose medium suitable for the growth of fastidious microorganisms.

For professional use only (in vitro diagnostic use).

Summary and Explanation

Traditionally blood agars have contained either casein hydrolysate or meat infusions. The advantage of the first lies in the production of large colonies and in the second in clearly defined haemolytic reactions and typical colony morphology. Ellner et al.¹ combined both types of peptone to create an improved medium, Columbia agar. Used without supplementation it is a good general-purpose isolation medium. Columbia agar is also used as a base for the addition of blood and in selective formulations when supplemented with antimicrobial agents. Columbia agar base with 5% sheep blood^{PLUS} has been developed for the cultivation of fastidious microorganisms.

Principle

Special peptone supplies nutrients for growth. Starch is added to absorb any toxic metabolites. The addition of sheep blood as well as acting as a nutrient allows the determination of the haemolytic properties of the isolate.

Typical Formula*

	grams per litre
Special peptone	23.0
Starch	1.0
Sodium chloride	5.0
Agar	10.0

Addition

	per litre
Defibrinated sheep blood	50.0 ml

* Adjusted as required to meet performance standards

Physical Characteristics

Colour	Traffic red
Clarity	Opaque
Fill weight	17 g ± 5 %
pH	7.3 ± 0.2

Precautions

This product is for *in vitro* diagnostic use and should only be used by trained individuals. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully.

For professional use only. Safety Data Sheet available on request.

Storage

This product is ready to use and no further preparation is necessary. Store product in its original packaging at 2–12°C until used.

Store away from light.

Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following the recommended guidelines.

Materials Required but Not Supplied

- (1) Inoculating loops, swabs, collection containers
- (2) Incubators
- (3) Quality control organisms

More information can be found at www.thermofisher.com

Procedure

Refer to local protocols and guidelines. The method followed will depend upon the sample under test. As an example of a protocol:

- (1) Inoculate the specimen onto the medium and streak to obtain isolated colonies.
- (2) If material is being cultured directly from a swab, roll the swab over a small surface area and streak for the isolation of single colonies.
- (3) Incubate the medium as appropriate for the sample type and microorganism under test.

The atmosphere, time and temperature of incubation will vary according to the type of specimen and the microorganism under test.

After incubation examine the plates for growth, colony morphology and haemolytic reactions. β haemolysis appears as a clear zone of surrounding the colony, α haemolysis appears as greenish pigmentation and a narrow zone of partial haemolysis surrounding the colony.

Quality Control

This medium can be tested with the following strains:

Incubation Conditions: 18 – 24 h @ 36 ± 1°C, aerobic

Positive Controls	
Inoculum of 50 – 120 colony forming units (cfu).	
<i>Streptococcus pyogenes</i> ATCC®12344	1 – 2 mm, light grey colonies with beta haemolysis
Colony count shall be greater than or equal to 50% of the control medium (Columbia agar with sheep blood).	
Inoculum of 10 ³ – 10 ⁴ cfu	
<i>Staphylococcus aureus</i> ATCC®25923	Good growth, white colonies
<i>Staphylococcus aureus</i> ATCC®6538	Good growth, yellow colonies with haemolysis
<i>Streptococcus pneumoniae</i> ATCC®6305	Good growth, dark grey colonies with alpha haemolysis
<i>Pseudomonas aeruginosa</i> ATCC®9027	Good growth, grey shiny colonies
<i>Escherichia coli</i> ATCC®8739	Good growth, dark grey colonies

Note:

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The product should not be used if:

- (1) The product is contaminated
- (2) The colour has changed
- (3) The expiration date has passed

(4) There are other signs of deterioration

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Limitations

Strains with specific growth requirements may not grow on this medium. Most strains of *Haemophilus* are unable to grow on sheep blood agars. Haemolytic reactions will vary by strain.

Packaging

PB5039A	Ten 90mm plates, film wrapped
----------------	-------------------------------

Bibliography

- (1) Ellner P.D., Stoessel C.J., Drakeford E. and Vasi F. (1966) *Tech. Bull. Reg. Med. Techn.* 36. No. 3, reprinted in *Amer. J. Clin. Path.* (1966) 45. 502-504.

Symbol Legend

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation (storage temp.)
	Use by (expiration date)
	Lot number
	Protect from light
	Consult instructions for use

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. ATCC® is a trademark of American Type Culture Collection. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

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Oxoid Deutschland GmbH