[www.thermofisher.com](http://www.thermofisher.com)

# Thermo

## SCIENTIFIC

Another common fungus is *Aspergillus brasiliensis*. *A. brasiliensis* is a species of mould and produces spores that can easily make their way into lungs and airways. It generally inhabits soil, and is so widespread indoors and outdoors that it is almost impossible to avoid breathing it in. This is especially dangerous for individuals with weakened immune systems, as it can quickly develop into lung and sinus infections. One of the most common infections is aspergillosis<sup>2</sup>. Aspergillosis is an opportunistic infection, and its severity can range from non-invasive and harmless to invasive and life-threatening. In immunocompromised patients, the infection tends to be highly invasive and fairly common. Invasive pulmonary aspergillosis is increasing in incidence among individuals with weakened immune systems, as well as in immunocompetent individuals, with high mortality for both<sup>2</sup>.

## Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol

**REF PO5096A, PO5243E\* and PO5258E\*\***

\* Biplate version with Chromogenic *C. albicans* Agar.  
\*\* Biplate version with *Brilliance™* Candida Agar.  
Additional IFUs available at [www.thermofisher.com](http://www.thermofisher.com).

### Intended Use **IVD**

Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol devices (PO5096A, PO5243E and PO5258E) are acidic pH selective media for the isolation of dermatophytes, other fungi, and yeasts from skin, hair, nails, genital, respiratory and urine samples from patients. Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having fungal infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

### Summary and Explanation

Dermatophytes are a specialised group of fungi that cause superficial fungal infections in the skin, hair and nails. Dermatophyte infections are arguably the most common form of fungal infection, with more than 20 different species. These species are divided into three genera: *Trichophyton*, *Microsporum* and *Epidermophyton*. Whilst the natural habitat of dermatophytes is soil, some these organisms have branched out and developed host specificity. Anthropophilic and zoophilic dermatophytes have the ability to infect humans. Common species of dermatophyte that cause fungal infections in humans are *Trichophyton rubrum*, *Trichophyton mentagrophytes* and *Microsporum canis*<sup>1</sup>. Yeasts are single-celled fungi that require warmth, moisture and nutrients to survive. Yeasts can survive normally inside the human body and on the skin without causing infections, however certain parts of the body have conditions suitable for yeast cultures to grow abundantly<sup>2</sup>.

Dermatophytes, yeasts and other fungi have the ability to cause invasive diseases, particularly in immunocompromised individuals. They are also the causative agents of superficial skin infections, such as those caused by *T. rubrum*, systemic infections, such as candidiasis, and chronic infections such as aspergillosis.

*Candida albicans* is a yeast that is a common commensal invader of the throat, skin, mouth, respiratory tract and genitals. Although *C. albicans* is a commensal organism in the human body, small changes in environment can lead to rapid growth, quickly resulting in infections such as candidiasis. This is especially common in immunocompromised individuals. One such example of this is when environmental changes occur in the vagina or the penis, causing rapid, uncontrolled proliferation of *C. albicans*. This results in vaginal candidosis or vaginitis in women and penile thrush in men<sup>3</sup>. Depending on the health of the individual, candidiasis can range from a superficial, mild fungal infection to systemic infections that have an increased risk in morbidity and mortality<sup>3,4-12</sup>. Although candidiasis can be caused by a number of *Candida* spp., the most common isolate from bloodstream infections, UTIs, burns and nail infections is *C. albicans*<sup>6, 9, 12</sup>.

Predisposition to fungal infections such as *Candida* and *Aspergillus* are caused by faults in the function of neutrophils, or lack of neutrophils. This is why the likelihood of an immunosuppressed individual becoming severely infected is much greater than individuals with normal white blood cell counts. As diagnosis is often delayed, due to failure in isolating or interpreting the presence of invasive fungi, or the presence of non-specific clinical symptoms, or the use of limited ability serotyping, fungal infections in immunocompromised hosts can be fatal<sup>13</sup>. Due to the severity of fungal infections in immunocompromised individuals, it is therefore highly important to isolate and identify dermatophytes, other fungi, and yeasts from skin, hair, nails, genital, respiratory and urine samples from immunocompromised patients. Early diagnosis is vital in the prevention and reduction of the morbidity and mortality associated with severe fungal infections. Presumptive identification of dermatophytes, yeasts and other fungi is carried out using a combination of microscopic appearance and cultural appearance on Sabouraud dextrose agar<sup>6</sup>.

### Principle of Method

Sabouraud Glucose Agar contains mycological peptone and glucose to support growth of fungi. It is slightly acidic which helps to inhibit bacteria whilst being favourable to fungi. Addition of antimicrobial agents, chloramphenicol and gentamicin increases the selectivity of the medium. Chloramphenicol is a bacteriostatic antimicrobial with a broad spectrum of activity against both Gram-positive and Gram-negative bacteria. Gentamicin is an aminoglycoside with activity against many Gram-negative aerobic bacteria and some strains of staphylococci.

### Typical Formula

	grams per litre
Mycological peptone	10.0
Glucose	40.0
Gentamicin	0.1
Chloramphenicol	0.05
Agar	5.0

### Physical Appearance

Colour	Ivory
Clarity	Transparent
Fill weight	17 ± 1.0g
pH	5.6 ± 0.2

### Materials Provided

PO5096A: 10 x 90 mm agar Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol plates.

PO5243E: 10 x 90 mm agar Sabouraud CG/Chromogenic *C. albicans* Agar biplates.

PO5258E: 10 x 90 mm *Brilliance™* Candida Agar/Sabouraud GC Agar biplates.

Each plate should only be used once.

## Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

## Storage

- Store product in its original packaging at 2°C to 12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

## Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product ([www.thermofisher.com](http://www.thermofisher.com)).

## Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

## Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 29 (Public Health England, 2020).

## Procedure

- Allow product to equilibrate to room temperature.
- Inoculate with the sample, for example skin dandruff, skin scrapings, hair and nail pieces, directly or after fragmentation and dilution.
- Incubate plates aerobically for 22 ± 1°C for 48 – 72 hours.
- Identify isolates by performing morphological, biochemical, and/or serological tests as required.

## Interpretation

White colonies (2 - 3 mm) indicate the presence of *Candida albicans*. *Aspergillus brasiliensis* presents as good growth, white mycelium, black spores.

## Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature, etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation conditions: 22 ± 1°C for 48 – 72 hours aerobic

Positive Control		
<i>Candida albicans</i> ATCC® 10231™	50 – 120 cfu*	2 - 3 mm white colonies
<i>Aspergillus brasiliensis</i> ATCC® 16404™	10 <sup>3</sup> – 10 <sup>4</sup> cfu** 1 colony***	Good growth, white mycelium, black spores
Negative Control		
<i>Escherichia coli</i> ATCC® 25922™ For PO5096A and PO5243E	≥ 10 <sup>4</sup> cfu	No growth
<i>Escherichia coli</i> ATCC® 25922™ For PO5258E	10 <sup>4</sup> -10 <sup>5</sup> cfu	No growth

\*Colony counts ≥ 50% of the control medium.

\*\*For PO5096A.

\*\*\*For PO5243E and PO5258E.

## Limitations

Some microorganisms, for example *Nocardia spp.* (filamentous bacteria) may be inhibited by chloramphenicol. Due to variation in nutritional requirements some strains may grow poorly or fail to grow, such as *Malassezia furfur* which requires a lipid supplement. Bacteria may grow if they are resistant to the antimicrobials present.

All identifications are presumptive and should be confirmed using appropriate methods. Medium choice and incubation protocol will depend upon the specimen type. A non-selective medium should be inoculated alongside the selective medium. Incubation at 30-37°C is suitable for yeasts, but dermatophytes may not grow at temperatures above 30°C.

## Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of dermatophyte, fungal and yeast strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices.

The precision of Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol devices has been demonstrated:

For Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol (PO5096A) - an overall pass rate of 100% over one month of testing (November – December 2021 – 10 batches).

For Sabouraud CG/Chromogenic *C. albicans* Agar (PO5243E) - an overall pass rate of 100% over six months of testing (June – December 2021 – 10 batches).

For *Brilliance Candida* Agar / Sabouraud GC Agar (PO5258E) - an overall pass rate of 100% over two months of testing (November 2021 – January 2022 – 10 batches).

The products have maintained a high pass rate, which helps demonstrate the precision of the devices, showing that the performance is reproducible.

The devices are tested in-house as part of the QC process since the products were launched in 1998 for Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol (PO5096A), 2004 for Sabouraud CG / Chromogenic *C. albicans* Agar (PO5243E) and 2012 for *Brilliance* Candida Agar / Sabouraud GC Agar (PO5258E).

For target organisms for Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol (PO5096A), when using 120 cfu and 10<sup>4</sup> cfu inoculum of *Candida albicans* and *Aspergillus brasiliensis* respectively and incubating the device at 21-23°C for 48-72 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria. For non-target organisms, when using ≥ 10<sup>4</sup> cfu inoculum of *Escherichia coli* and incubating the device at 21-23°C for 48-72 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria.

For target organisms for Sabouraud CG / Chromogenic *C. albicans* Agar (PO5243E), on the Sabouraud GC Agar side, when using 10<sup>4</sup> cfu inoculum of *Candida albicans* and the 1 colony direct streaking method for *Aspergillus brasiliensis* respectively and incubating the device at 21-23°C for 48-72 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria. On the Chromogenic *C. albicans* Agar side, when using 10<sup>4</sup> cfu inoculum of both *Candida albicans* and *Candida tropicalis* and incubating the device at 21-23°C for 48-72 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria. For non-target organisms, when using ≥ 10<sup>4</sup> cfu inoculum of *Escherichia coli* on both sides of the biplate and incubating the device at 21-23°C for 48-72 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria.


For target organisms for *Brilliance* Candida Agar / Sabouraud GC Agar (PO5258E), on the Sabouraud GC Agar side, when using 10<sup>4</sup> cfu inoculum of *Candida albicans* and the 1 colony direct streaking method for *Aspergillus brasiliensis* respectively and incubating the device at 31-33°C for 48-72 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria. On the *Brilliance* Candida Agar side, when using 10<sup>4</sup> cfu inoculum of both *Candida albicans* and *Candida tropicalis* and incubating the device at 31-33°C for 48-72 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria. For non-target organisms, when using 10<sup>5</sup> cfu inoculum of *Escherichia coli* on both sides of the biplate and incubating the device at 31-33°C for 48-72 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria.

### Bibliography

1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'.

3. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
4. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
5. Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
6. Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
7. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
8. Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
9. Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
10. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
11. Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2). <https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
12. Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
13. Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
14. Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. *Clinical Haematology* 13 (3): 599-630. <https://pubmed.ncbi.nlm.nih.gov/6388935/>

**Symbol Legend**

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limitation (storage temp.)
	Use by (expiration date) YYYY-MM
	Keep away from sunlight
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Manufacturer
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Assessment
	UK Conformity Assessment
Made in Germany	Made in Germany

© 2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Oxoid GmbH, Am Lippeglacis 4-8, 46483  
Wesel, Germany

For technical assistance please contact your local distributor.

**Revision information**

Version	Date of issue and modifications introduced
1.0	2022-08-04. New document. (LIVE)

[www.thermofisher.com](http://www.thermofisher.com)

**Thermo**  
SCIENTIFIC

## Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol

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## Bibliographie

1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
3. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
4. Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
5. Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
6. Public Health England. 2016. Investigation of Dermatological Specimens for Superficial

- Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1).  
<https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2).  
<https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
  - Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1).  
<https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
  - Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5).  
<https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
  - Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2).  
<https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
  - Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5).  
<https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
  - Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7).  
<https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
  - Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599–630.  
<https://pubmed.ncbi.nlm.nih.gov/6388935/>

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### Bibliographie

1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
3. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
4. Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
5. Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
6. Public Health England. 2016. Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.

7. Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2). <https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
11. Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
12. Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
13. Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599–630. <https://pubmed.ncbi.nlm.nih.gov/6388935/>

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1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
3. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6).

- <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1).  
<https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
  - Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9).  
<https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
  - Public Health England. 2016. Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1).  
<https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
  - Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2).  
<https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
  - Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1).  
<https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
  - Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5).  
<https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
  - Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2).  
<https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
  - Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5).  
<https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
  - Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7).  
<https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
  - Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599-630.  
<https://pubmed.ncbi.nlm.nih.gov/6388935/>

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## Bibliografia

1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH\\_86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH_86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
3. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
4. Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
5. Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
6. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
7. Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5).

- <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2).  
<https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
11. Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5).  
<https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
12. Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7).  
<https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
13. Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599–630.  
<https://pubmed.ncbi.nlm.nih.gov/6388935/>

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## Bibliografia

1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
3. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
4. Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
5. Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
6. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
7. Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2).

- <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1).  
<https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
  - Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5).  
<https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
  - Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2).  
<https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
  - Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5).  
<https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
  - Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7).  
<https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
  - Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599–630.  
<https://pubmed.ncbi.nlm.nih.gov/6388935/>

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## Bibliografia

1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
3. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
4. Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
5. Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
6. Public Health England. 2016. Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
7. Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.

10. Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2). <https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
11. Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
12. Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
13. Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599–630. <https://pubmed.ncbi.nlm.nih.gov/6388935/>

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## Bibliografie

1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
3. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
4. Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
5. Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
6. Public Health England. 2016. Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
7. Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than

- Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2). <https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
11. Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
12. Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
13. Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599–630. <https://pubmed.ncbi.nlm.nih.gov/6388935/>

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## Zdroje

- World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Accessed online January 4 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
- Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
- Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2). <https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
- Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for

- Microbiology Investigations B 57 (3.5).  
<https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
12. Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7).  
<https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
13. Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599–630.  
<https://pubmed.ncbi.nlm.nih.gov/6388935/>

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1. World Health Organisation. 1987. 'Guidelines for the Diagnosis, Prevention and Control of Dermatophytoses in Man and Animals'. Consultado en línea el 4 de enero de 2022. [https://apps.who.int/iris/bitstream/handle/10665/61519/WHO\\_CDS\\_VPH.86.67.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/61519/WHO_CDS_VPH.86.67.pdf?sequence=1&isAllowed=y).
2. Centers for Disease Control and Prevention. 2021. 'Vaginal Candidiasis'. <https://www.cdc.gov/fungal/diseases/candidiasis/genital/index.html>.  
———. 2021. 'Where Aspergillosis Comes From'. <https://www.cdc.gov/fungal/diseases/aspergillosis/causes.html>.  
———. 2021. 'Information for Healthcare Professionals about Aspergillosis'. <https://www.cdc.gov/fungal/diseases/aspergillosis/health-professionals.html>.
3. Public Health England. 2017. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
4. Public Health England. 2015. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
5. Public Health England. 2015. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
6. Public Health England. 2016. Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
7. Public Health England. 2017. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi>

- id-1-preliminary-identification-of-medically-important-bacteria.
- Public Health England. 2017. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1).  
<https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
  - Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5).  
<https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
  - Public Health England. 2019. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2).  
<https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
  - Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5).  
<https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
  - Public Health England. 2019. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7).  
<https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
  - Hawkins, C, D Armstrong. 1984. 'Fungal Infections in the Immunocompromised Host'. Clinical Haematology 13 (3): 599–630.  
<https://pubmed.ncbi.nlm.nih.gov/6388935/>

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Oxoid GmbH, Am Lippeglacis 4-8, 46483 Wesel, Alemania