## Oxoid Prepared Medium

### Oxoid Staph / Strep Selective Medium (Columbia CNA)

REF PB5049A

#### Intended Use

IVD

ThermoScientific™ Oxoid™ Staph / Strep Selective Medium (Columbia CNA) is a selective medium for the isolation of staphylococci and streptococci.

For professional use only (in vitro diagnostic use).

#### Summary and Explanation

Traditionally blood agars have contained either casein hydrolysate or meat infusions. The advantage of the first lies in the production of large colonies and in the second in clearly defined haemolytic reactions and typical colony morphology. Ellner et al.<sup>1</sup> combined both types of peptone to create Columbia agar an improved all-round medium. Used without supplementation it is good general-purpose medium. Columbia agar is also used as a base medium for the addition of blood and in selective formulations when supplemented with antimicrobial agents. Staph / Strep Selective medium is used for the isolation of streptococci and staphylococci. Colistin and nalidixic acid were suggested as selective agents for the isolation of Gram-positive cocci by Ellner et al.<sup>1</sup>. Important Gram-positive cocci are more readily recognized and easily isolated from the mixed bacterial populations contained in many food and clinical samples.

#### Principle

Special peptone supplies nutrients for growth. Starch is added to absorb any toxic metabolites. Sheep blood is a nutrient source and allows the determination of the haemolytic properties of the isolate. Colistin sulphate is a polymixin antibiotic that is active against Gram-negative bacteria including *Pseudomonas aeruginosa*; *Proteus* spp., however, are not inhibited. Nalidixic acid is also active mainly against Gram-negative bacteria including *Proteus* spp. The Oxoid medium is modified from that of Ellner<sup>1</sup> in that it contains a reduced concentration of the antimicrobial agents.

#### Typical Formula\*

	grams per litre
Special peptone	23.0
Starch	1.0
Sodium chloride	5.0
Agar	10.0

#### Additions

	per litre
Defibrinated sheep blood	70.0 mL
Nalidixic acid	0.005 g
Colistin	0.0075 g

\* Adjusted as required to meet performance standards

#### Physical Characteristics

Colour	Traffic red
Clarity	Opaque
Fill weight	17 g ± 5 %
pH	7.3 ± 0.2

#### Precautions

This product is for *in vitro* diagnostic use and should only be used by trained individuals. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products. It

is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully.

For professional use only. Safety Data Sheet available on request.

#### Storage

This product is ready to use and no further preparation is necessary. Store product in its original packaging at 2–12°C until used. Store away from light. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

#### Specimen Collection, Handling and Storage

Specimens should be collected and handled following the recommended guidelines.

#### Materials Required but Not Supplied

- (1) Inoculating loops, swabs, collection containers
- (2) Incubators
- (3) Quality control organisms

More information can be found at [www.thermofisher.com](http://www.thermofisher.com)

#### Procedure

Refer to any local protocols and guidelines.

- (1) Inoculate the specimen onto the medium and streak to obtain isolated colonies.
- (2) If material is being cultured directly from a swab, roll the swab over a small surface area and streak for the isolation of single colonies.
- (3) Incubate the medium at 36 ± 1°C for 24 hours to 48 hours aerobically.

After incubation examine the plates for growth, colony morphology and haemolytic reactions.

#### Quality Control

This medium can be tested with the following strains:

Incubation Conditions: 24 h @ 36 ± 1°C, aerobic

Positive Controls	
Inoculum of 50 - 120 colony forming units (cfu).	
<i>Streptococcus pyogenes</i> ATCC®12344	1 – 2 mm, grey shiny colonies with β-haemolysis
Colony count shall be greater than or equal to 50% of the control medium (Columbia agar with sheep blood).	
Inoculum of 10 <sup>3</sup> - 10 <sup>4</sup> cfu	
<i>Staphylococcus aureus</i> ATCC®25923	Good growth White shiny colonies
<i>Streptococcus pneumoniae</i> ATCC®6305	Good growth, grey colonies with α-haemolysis
Negative Controls	
Inoculum of ≥ 10 <sup>4</sup> cfu	
<i>Escherichia coli</i> ATCC®25922	Complete inhibition (≤ 10 cfu)
<i>Pseudomonas aeruginosa</i> ATCC®27853	Complete inhibition (≤ 10 cfu)

#### Note:

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The product should not be used if:

- (1) The product is contaminated
- (2) The colour has changed
- (3) The expiration date has passed
- (4) There are other signs of deterioration

**Performance**

Performance was evaluated using 30 bacterial strains including the following; *Staphylococcus aureus*, *S. pneumoniae*, *Streptococcus viridans*, *S. pyogenes*, *Streptococcus agalactiae*, *E.coli*, *Pseudomonas aeruginosa*, *Proteus mirabilis*, *Proteus hauseri*, *Enterococcus faecalis*, *Enterococcus faecium*, *Moraxella catarrhalis*, *Clostridium perfringens*, *Clostridium sporogenes*, *Pasteurella multocida*, and *Bacteroides fragilis*. All organisms gave expected growth characteristics according to the current product specification.

**Limitations**

The haemolytic reaction will be strain dependent; it may also vary for some microorganisms with the type of blood used. All identifications are presumptive and should be confirmed using appropriate methods. Due to variation in nutritional requirements or sensitivity to selective agents some strains of the target organisms may be encountered that grow poorly or fail to grow on this medium. Incubation in a carbon dioxide enriched atmosphere may cause inhibition of staphylococci<sup>2</sup>. Microorganisms other than streptococci and staphylococci which are resistant to the antimicrobial agents used in the medium will be able to grow.

**Packaging**

<b>PB5049A</b>	Ten 90mm plates, film wrapped
----------------	-------------------------------

**Bibliography**

- (1) Ellner P.D., Stoessel C.J., Drakeford E. and Vasi F. (1966) *Tech. Bull. Reg. Med. Techn.* 36. No. 3, reprinted in *Amer. J. Clin. Path.* (1966) 45. 502-504.
- (2) Morton C.E.G., Holt H.A. (1989) *Medical Laboratory Sciences*. 46(1): 72-73.

**Symbol Legend**

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation (storage temp.)
	Use by (expiration date)
	Lot number
	Protect from light
	Consult instructions for use

trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. ATCC® is a trademark of American Type Culture Collection. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

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Version 4

 The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalog marks are

## Oxoid Fertignährmedium

### Oxoid Staph / Strep Selective Medium (Columbia CNA)

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Oxoid Deutschland GmbH