www.thermofisher.com

## Mueller Hinton Agar with Horse Blood

**REF** PB5303A

### Intended Use

Mueller-Hinton Agar with Horse Blood is an antimicrobial susceptibility agar recommended for disc diffusion and MIC testing against fastidious microorganisms isolated from clinical samples.

The device is for professional use only, is not automated nor is it a companion diagnostic.

### Summary and Explanation

Streptococci are Gram-positive cocci which inhabit a wide range of ecological niches, from non-pathogenic species used in the yoghurt making industry, to commensal members of the mammal oral microbiome, and invasive human pathogens<sup>1,2,3</sup>. Streptococci can be classified by serotyping based on cell-wall carbohydrate composition into the Lancefield groups A to S, and the non-Lancefield Streptococci<sup>1,2,3</sup>. Non-Lancefield Streptococci cannot be grouped by Lancefield serotyping, and include a broad range of non-pathogenic and pathogenic species with *S. pneumoniae* thought to be the most significant human pathogen<sup>1,2,3</sup>.

A number of National Institute for Health and Care Excellence (NICE) and Public Health England (PHE) guidelines describe *S. pneumoniae* as one of the leading causes of meningitis and bloodstream infections, both invasive diseases associated with high morbidity and mortality<sup>3,4,5,6,7,8</sup>. *S. pneumoniae* is a notifiable organism, meaning that infection with this pathogen is of such clinical significance that it is required by law to report diagnosis of these diseases to the government authority PHE<sup>3,9</sup>. The status of *S. pneumoniae* as a notifiable organism highlights the importance of infection surveillance to identify future disease outbreaks and reduce potential morbidity and mortality<sup>10,11</sup>.

### Principle of Method

Mueller-Hinton Agar has been criticised because of variation in performance between manufacturers and batches of medium. These effects were shown as:

- (1) Minimum inhibitory concentration (MIC) variations with aminoglycosides against *Pseudomonas aeruginosa* and tetracycline against *Staphylococci*. This may have been due to differing concentrations of the divalent cations of calcium and magnesium.
- (2) Variation in thymine and thymidine content, which affect sulphonamide and trimethoprim MIC values
- (3) Differences in the characteristics of the agar used in the medium, especially diffusion properties

In light of the criticisms the Clinical Laboratory Standards Institute (CLSI) (previously NCCLS) called interested manufacturers together to discuss the standardisation and stabilisation of Mueller-Hinton Agar. Control methods were established whereby critical antimicrobial/organism combinations had to yield consistent zones of inhibition within two mm of the specified diameters in the standard. The result of this cooperative effort is that Mueller-Hinton Agar is now a standard medium.

# Thermo

## SCIENTIFIC

Supplemented Mueller-Hinton Agar is recommended both by CLSI and European Committee on Antimicrobial Susceptibility Testing (EUCAST) for testing non-fastidious strains. Supplements are added to make it suitable for more fastidious strains.

### Typical Formula

	grams per litre
Beef, dehydrated infusion from	300.0
Casein hydrolysate	17.5
Starch	1.5
NAD	0.02
Agar	17.0
Defibrinated horse blood	50.0ml

### Physical Appearance

Colour	Signal red
Clarity	Opaque
Fill weight	23.5 ± 5%
pH	7.3 ± 0.1

### Materials Provided

PB5303A: 10 x 90mm Mueller Hinton agar with horse blood plates

Each plate should only be used once.

### Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

### Optional Materials

- Erythromycin antimicrobial susceptibility discs (CT0020B)
- Cefpodoxime antimicrobial susceptibility discs (CT1612B)
- Ertapenem antimicrobial susceptibility discs (CT1761B)
- Moxifloxacin antimicrobial susceptibility discs (CT1633B)

### Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

### Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused

reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products. Refer to the Material Safety Data Sheet (MSDS) for safe handling and disposal of the product ([www.thermofisher.com](http://www.thermofisher.com)).

### Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

### Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the CLSI standard M100<sup>10</sup>.

### Procedure

Allow plates and antimicrobial susceptibility disks to equilibrate to room temperature before use. The agar surface should not have excess moisture prior to inoculation.

Inoculum may be prepared utilizing either the growth method or the direct colony method. Refer to appropriate DIN/CLSI documents for guidelines specific to each test organism

#### CLSI Method:

1. Growth method:
  - a. Select at least 3-5 isolated colonies of the same morphology type from the agar culture. Touch the top of each colony with a loop and transfer to 4-5ml of a suitable broth medium such as Tryptic Soy Broth (TSB)
  - b. Incubate the broth culture at 35-37°C until the turbidity equals or exceeds that of a 0.5 McFarland or equivalent.
  - c. If necessary, adjust the turbidity of the suspension with sterile saline or broth to achieve a turbidity to a 0.5 McFarland standard. This step may be performed visually or using a photometric device. If performed visually use adequate light to compare the suspension to the McFarland 0.5 against a card with a white background and contrasting black lines<sup>5</sup>.
2. Direct Colony Method: (Method of choice for streptococci and staphylococci).
  - a. Prepare a direct suspension of the test organism in saline or broth from an 18-24 hour culture on nonselective media.
  - b. Adjust the turbidity of the suspension as described under Growth method.
3. Inoculate Agar plates within 15 minutes of preparing organism suspension.
4. Immerse a sterile swab into the suspension. Rotate the swab against the side of the tube above the fluid level to remove excess fluid.
5. Inoculate the surface in three planes by rotating the plate approximately 60 degrees each time.
6. Replace the lid and allow the plate to rest on the bench at least 3 minutes but no longer than 15 minutes for the inoculum to be absorbed before applying antimicrobial susceptibility disks.
7. Apply disks individually or by using an antimicrobial disk dispenser. Disk should not be closer than 24 mm from centre to centre and not placed too close to the edge of the plate. Because some drug diffuses almost instantaneously, do not relocate a disk once it has come in contact with the agar surface. Tap disks gently with sterile needle or forceps to ensure complete contact with the agar surface.

8. Invert the plate and place in the incubator within 15 minutes of disk application.
  - a. Incubate Mueller-Hinton Agar with Sheep Blood plates aerobically at 35-37°C for 16-18 hours.
 

**Note:** Staphylococci and enterococci require a full 24 hours incubation. Methicillin-resistant staphylococci (MRS) may not be detected at incubation temperatures above 35°C.
  - b. Incubate *Streptococcus* spp. in 5-7% CO<sub>2</sub> at 35-37°C for 20-24 hours.
  - c. Refer to CLSI document M45 for incubation times, temperatures, and atmospheres for fastidious or infrequently isolated bacteria<sup>19</sup>.

### Interpretation

After 16-24 hours incubation at 35–37°C, *S. pneumoniae* typically exhibit 1-2 mm colonies and exhibit α-haemolysis on blood agar. *S. pneumoniae* can also appear as draughtsman colonies with a depressed centre and elevated rim.

### Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

Incubation Conditions: 18 ± 2h @ 36 ± 1°C in an enhanced carbon dioxide atmosphere

Microorganism	Antibiotic disc	Zone of Inhibition Diameter (mm)
<i>Streptococcus pneumoniae</i> ATCC®49619	Erythromycin (E15)	26-32
	Cefpodoxime (CPD 10)	29-35
	Ertapenem (ETP 10)	28-34
	Moxifloxacin (MXF 5)	24-30
<i>Haemophilus influenzae</i> ATCC®49766	Erythromycin (E15)	10-16
	Cefpodoxime (CPD 10)	30-36
	Moxifloxacin (MXF 5)	30-36

### Limitations

Organisms with atypical enzyme patterns may give anomalous reactions on Mueller-Hinton Agar with Horse Blood.

### Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of fastidious microorganisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of Mueller-Hinton Agar with Horse Blood (PB5303A) was demonstrated by an overall pass rate of 99.7% obtained for the product over an 11 year period (2010-2021; 3439 batches) This shows that the performance is reproducible.

Performance is assessed as the zone of inhibition diameter measured on the test medium using Thermo Scientific™ Oxoid™ antimicrobial susceptibility discs, compared with the specified diameters in the CLSI standard M100<sup>10</sup>.

The inoculum is prepared from working cultures, diluted to obtain the specified inoculum levels and the test media is inoculated. Antibiotic discs are loaded in the same order and concentrations as they appear in the inspection plan for

the medium. Following incubation at 34-36°C for 16-20 hours in 4-6% CO<sub>2</sub>, the zone of inhibition diameter of the test organisms is measured, and the results and conclusions can be created.

**Bibliography**

1. Patterson, Maria Jevitz. 1996. "Streptococcus." In Medical Microbiology. 4th Edition. Galveston (TX): University of Texas Medical Branch at Galveston.
2. Dion, Christopher F., and John V. Ashurst. 2020. "Streptococcus Pneumoniae." In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
3. Public Health England. 2014a. "Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms." UK Standards for Microbiology Investigations.
4. Public Health England. 2014b. "Investigation of Ear Infections and Associated Specimens." UK Standards for Microbiology Investigations.
5. National Institute for Health and Care Excellence. 2019b. "Fever in under 5s: Assessment and Initial Management."
6. National Institute for Health and Care Excellence. 2010. "Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management."
7. National Institute for Health and Care Excellence. 2016. "SepsiTest Assay for Rapidly Identifying Bloodstream Bacteria and Fungi."
8. Public Health England. 2016. "Investigation of Pus and Exudates." UK Standards for Microbiology Investigations.
9. Public Health England. 2019. "Statutory Notification by Registered Medical Practitioners of All Hazards: Infections, Chemicals & Radiation."
10. Clinical & Laboratory Standards Institute. 2020. "Performance Standards for Antimicrobial Susceptibility Testing." <https://clsi.org/standards/products/microbiology/documents/m100/>.

	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. CLSI is a trademark of the Clinical Laboratory and Standards Institute. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

**Symbol Legend**

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use

 Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, England

For technical assistance please contact your local distributor.

**Revision information**

Version	Date of issue and modifications introduced
1.0	2022-11-02. New document. (LIVE)

www.thermofisher.com

## Mueller Hinton Agar with Horse Blood

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## Bibliographie

1. Patterson, Maria Jevitz. 1996. "Streptococcus." In Medical Microbiology. 4th Edition. Galveston (TX): University of Texas Medical Branch at Galveston.
2. Dion, Christopher F., and John V. Ashurst. 2020. "Streptococcus Pneumoniae." In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
3. Public Health England. 2014a. "Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms." UK Standards for Microbiology Investigations.
4. Public Health England. 2014b. "Investigation of Ear Infections and Associated Specimens." UK Standards for Microbiology Investigations.
5. National Institute for Health and Care Excellence. 2019b. "Fever in under 5s: Assessment and Initial Management."
6. National Institute for Health and Care Excellence. 2010. "Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management."
7. National Institute for Health and Care Excellence. 2016. "SepsisTest Assay for Rapidly Identifying Bloodstream Bacteria and Fungi."
8. Public Health England. 2016. "Investigation of Pus and Exudates." UK Standards for Microbiology Investigations.
9. Public Health England. 2019. "Statutory Notification by Registered Medical Practitioners of All Hazards: Infections, Chemicals & Radiation."
10. Clinical & Laboratory Standards Institute. 2020. "Performance Standards for Antimicrobial Susceptibility Testing." <https://clsi.org/standards/products/microbiology/documents/m100/>.

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 Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Angleterre

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### Bibliographie

1. Patterson, Maria Jevitz. 1996. "Streptococcus." In Medical Microbiology. 4th Edition. Galveston (TX): University of Texas Medical Branch at Galveston.
2. Dion, Christopher F., and John V. Ashurst. 2020. "Streptococcus Pneumoniae." In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
3. Public Health England. 2014a. "Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms." UK Standards for Microbiology Investigations.
4. Public Health England. 2014b. "Investigation of Ear Infections and Associated Specimens." UK Standards for Microbiology Investigations.
5. National Institute for Health and Care Excellence. 2019b. "Fever in under 5s: Assessment and Initial Management."
6. National Institute for Health and Care Excellence. 2010. "Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management."
7. National Institute for Health and Care Excellence. 2016. "SepsisTest Assay for Rapidly Identifying Bloodstream Bacteria and Fungi."
8. Public Health England. 2016. "Investigation of Pus and Exudates." UK Standards for Microbiology Investigations.

9. Public Health England. 2019. "Statutory Notification by Registered Medical Practitioners of All Hazards: Infections, Chemicals & Radiation."
10. Clinical & Laboratory Standards Institute. 2020. "Performance Standards for Antimicrobial Susceptibility Testing." <https://clsi.org/standards/products/microbiology/documents/m100/>.

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Oxoid Limited, Wade Road, Basingstoke,  
Hampshire, RG24 8PW, England

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1. Patterson, Maria Jevitz. 1996. “Streptococcus.” In Medical Microbiology. 4th Edition. Galveston (TX): University of Texas Medical Branch at Galveston.
2. Dion, Christopher F., and John V. Ashurst. 2020. “Streptococcus Pneumoniae.” In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
3. Public Health England. 2014a. “Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms.” UK Standards for Microbiology Investigations.
4. Public Health England. 2014b. “Investigation of Ear Infections and Associated Specimens.” UK Standards for Microbiology Investigations.
5. National Institute for Health and Care Excellence. 2019b. “Fever in under 5s: Assessment and Initial Management.”
6. National Institute for Health and Care Excellence. 2010. “Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management.”
7. National Institute for Health and Care Excellence. 2016. “SepsisTest Assay for Rapidly Identifying Bloodstream Bacteria and Fungi.”
8. Public Health England. 2016. “Investigation of Pus and Exudates.” UK Standards for Microbiology Investigations.
9. Public Health England. 2019. “Statutory Notification by Registered Medical Practitioners

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Oxoid Limited, Wade Road, Basingstoke,  
Hampshire, RG24 8PW, H.B.

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- National Institute for Health and Care Excellence. 2016. "SepsiTest Assay for Rapidly Identifying Bloodstream Bacteria and Fungi."
- Public Health England. 2016. "Investigation of Pus and Exudates." UK Standards for Microbiology Investigations.
- Public Health England. 2019. "Statutory Notification by Registered Medical Practitioners of All Hazards: Infections, Chemicals & Radiation."
- Clinical & Laboratory Standards Institute. 2020. "Performance Standards for Antimicrobial Susceptibility Testing." <https://clsi.org/standards/products/microbiology/documents/m100/>.

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### Bibliografia

- Patterson, Maria Jevitz. 1996. "Streptococcus." In Medical Microbiology. 4th Edition. Galveston (TX): University of Texas Medical Branch at Galveston.
- Dion, Christopher F., and John V. Ashurst. 2020. "Streptococcus Pneumoniae." In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
- Public Health England. 2014a. "Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms." UK Standards for Microbiology Investigations.
- Public Health England. 2014b. "Investigation of Ear Infections and Associated Specimens." UK Standards for Microbiology Investigations.
- National Institute for Health and Care Excellence. 2019b. "Fever in under 5s: Assessment and Initial Management."
- National Institute for Health and Care Excellence. 2010. "Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management."

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Oxoid Limited, Wade Road, Basingstoke,  
Hampshire, RG24 8PW, Anglia

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1. Patterson, Maria Jevitz. 1996. "Streptococcus." In Medical Microbiology. 4th Edition. Galveston (TX): University of Texas Medical Branch at Galveston.
2. Dion, Christopher F., and John V. Ashurst. 2020. "Streptococcus Pneumoniae." In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
3. Public Health England. 2014a. "Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms." UK Standards for Microbiology Investigations.
4. Public Health England. 2014b. "Investigation of Ear Infections and Associated Specimens." UK Standards for Microbiology Investigations.
5. National Institute for Health and Care Excellence. 2019b. "Fever in under 5s: Assessment and Initial Management." UK Standards for Microbiology Investigations.
6. National Institute for Health and Care Excellence. 2010. "Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management." UK Standards for Microbiology Investigations.
7. National Institute for Health and Care Excellence. 2016. "SepsiTest Assay for Rapidly Identifying Bloodstream Bacteria and Fungi." UK Standards for Microbiology Investigations.
8. Public Health England. 2016. "Investigation of Pus and Exudates." UK Standards for Microbiology Investigations.
9. Public Health England. 2019. "Statutory Notification by Registered Medical Practitioners of All Hazards: Infections, Chemicals & Radiation."
10. Clinical & Laboratory Standards Institute. 2020. "Performance Standards for Antimicrobial Susceptibility Testing." <https://clsi.org/standards/products/microbiology/documents/m100/>.

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