Mueller Hinton Agar

[REF] PO5007A, PO5403I

Intended Use

Mueller Hinton Agar (PO5007A and PO5403I) devices are antimicrobial susceptibility agar devices recommended for disc diffusion and Minimum Inhibitory Concentration (MIC) testing against isolated fastidious microorganisms.

The medium has been developed according to EUCAST (European Committee on Antimicrobial Susceptibility Testing) and CLSI (Clinical Laboratory and Standards Institute, Inc.) recommendations. Mueller Hinton Agar (PO5007A and PO5403I) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated and nor are they companion diagnostics.

Summary and Explanation

The vast majority of infectious diseases are caused by bacteria. The development of laboratory methods to grow these microorganisms using appropriate growth media has also allowed the determination of the sensitivity and resistance of pathogens to a variety of antimicrobial agents, allowing healthcare professionals to implement appropriate treatment regimens. The rapid global emergence of resistant bacteria threatens the effectiveness of antibiotics, which have revolutionised medicine and saved millions of lives. A number of bacteria have been labelled as showing urgent, serious, and worrying threats¹.

Assays for determination of antimicrobial susceptibility can be performed using genotypic or phenotypic techniques. The MIC serves as the 'gold standard' phenotypic resistance determination technique, with the organism classified as susceptible, intermediate, or resistant to the antimicrobial agent based on clinical MIC breakpoints. MIC assays, however, can be labour-intensive, time-consuming and require experience to interpret. Disc diffusion is widely used in routine clinical microbiology laboratories as it is low cost, simple to perform, versatile and almost all antimicrobial agents and pathogens can be tested. Disc diffusion is accurate and reproducible when carried out in strict accordance with guidelines/recommendations (such as EUCAST) or CLSI).

Principle of Method

Mueller Hinton Agar (PO5007A and PO5403I) devices contain dehydrated infusion from beef and casein hydrolysate to provide nitrogen, vitamins, carbon, and amino acids. Soluble starch is added to absorb toxic metabolites produced. Agar is added as the solidifying agent.

Mueller Hinton Agar devices are stable, reproducible, and standardized antimicrobial susceptibility agar devices recommended for disc diffusion and MIC testing.

Mueller Hinton Agar devices may be supplemented and adjusted to conform to Clinical and Laboratory Standards Institute (CLSI) and European Committee on Antimicrobial Susceptibility Testing (EUCAST) performance standards. Mueller Hinton Agar devices will require supplementation with the divalent cations Mg⁺⁺ and Ca⁺⁺ after sterilisation. The CLSI and EUCAST recommend the following cation levels Ca⁺⁺ 20-25mg/litre; Mg⁺⁺ 10-12.5mg/litre. For testing daptomycin, CLSI and EUCAST also recommend raising calcium to 50 mg/L. The medium needs to have a

zinc concentration lower than 3 mg/l to prevent resistant interpretations when testing carbapenems.

Typical Formula

Typical Formulation*	grams per litre (g/L)
Beef, dehydrated infusion from	300.0
Casein hydrolysate	17.5
Starch	1.5
Agar	17.0

Physical Appearance

Colour	Ivory
Clarity	Transparent
Filling Weight (PO5007A)	23.5 ± 5%
Filling Weight (PO5403I)	60 ± 1.2g
pH	7.3 ± 0.1

Materials Provided

PO5007A: 10 x 90mm Mueller Hinton agar plates
PO5403I: 4 x 120mm Mueller Hinton agar plates

Each plate should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubator
- Quality control organisms
- Antimicrobial susceptibility discs

Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the

relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Microorganisms isolated from clinical samples are tested on Mueller Hinton Agar (PO5007A and PO5403I) devices. Clinical samples are not inoculated directly onto the medium

Procedure

- Allow product to equilibrate to room temperature.
- Perform testing in strict accordance with current guidelines/recommendations in order to maintain acceptable levels of accuracy.

Interpretation

Antibiotic susceptibility tests should be performed and interpreted in accordance with current guidelines (such as those published by EUCAST, CLSI or ISO).

Quality Control

This medium can be tested with the following strains:
Incubation Conditions: 16–24 hours aerobically at 36±1°C
Incubation conditions: 18 ± 2 h @ 36 ± 1°C, aerobic

<i>Escherichia coli</i> ATCC® 25922™	
Antimicrobial	Zone Diameter (mm)
Amoxicillin/Clavulanic Acid (AMC 30)	18-24
Ciprofloxacin (CIP 5)	29-37
Gentamicin (CN 10)	19-36
Trimethoprim/ Sulphamethoxazole (SXT 25)	23-29

<i>Enterococcus faecalis</i> ATCC® 29212™	
Antimicrobial	Zone Diameter (mm)
Nitrofurantoin (F 100)	18-24
Ampicillin (AMP 2)	15-21
Ciprofloxacin (CIP 5)	19-25

<i>Staphylococcus aureus</i> ATCC® 29213™	
Antimicrobial	Zone Diameter (mm)
Gentamicin (CN10)	19-25
Tetracycline (TE30)	23-31
Penacillin G (P1)	12-18
Cefoxitin (FOX30)	24-30

<i>Pseudomonas aeruginosa</i> ATCC® 27853™	
Antimicrobial	Zone Diameter (mm)
Aztreonam (ATM30)	23-29
Imipenem (IPM10)	20-28
Cefepime (FEP30)	25-31

Limitations

The media may not support the growth of some fastidious strains, for example microorganisms that require thymine and thymidine for growth. Close adherence to the methodology is required to obtain satisfactory results.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct MIC or disc diffusion testing results are confirmed by the inclusion of well-characterised isolates and antimicrobials in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Mueller Hinton Agar (square plates)

(PO5403I) was demonstrated by an overall pass rate of 100 % obtained for the product over nine months of testing (12.10.2021 – 25.07.2022; 10 batches). The precision of Mueller Hinton Agar (PO5007A) was demonstrated by an overall pass rate of 100 % obtained for the product over one month of testing (25.07.2022 – 18.08.2022; 10 batches).


Mueller Hinton Agar (square plates) (PO5403I) and Mueller Hinton Agar (PO5007A) have been tested in-house as part of the QC process since they were first manufactured in 2003. When using optimum microorganism inoculum levels and incubating at the recommended time and temperature conditions, the user can achieve zone size diameters within the range stated in the specification, and inline with the EUCAST and/or CLSI guidelines.

Bibliography

1) Centers for Disease Control and Prevention (CDC). 2013. 'Office of Infectious Disease. Antibiotic resistance threats in the United States'. <http://www.cdc.gov/drugresistance/threat-report-2013>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment

	UK Conformity Assessment
	Unique device identifier
	Made in Germany

©2022 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
 CLSI is a trademark of the Clinical Laboratory and Standards Institute.
 All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

Version	Date of modifications introduced
1.0	2022-10-22 Original Document

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©2022 Thermo Fisher Scientific Inc. Todos los derechos reservados.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus filiales.
Oxoid GmbH Deutschland,